
Suprasorb® CNP endo Kit Oral

Zur Unterdrucktherapie im Ösophagus.
Pour traitement par pression négative dans l'œsophage.
For negative pressure therapy in the oesophagus.
Per la terapia a pressione negativa nell'esofago.

Assembled by:
Lohmann & Rauscher International GmbH & Co. KG
Westerwaldstraße 4 · 56579 Rengsdorf, Germany

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en Instructions for Use

Table of Contents

- 1) Kit description and composition
- 2) Required material
- 3) Intended purpose
- 4) Indications
- 5) Contraindications
- 6) Side effects
- 7) Interactions
- 8) Storage
- 9) Preparation
- 10) Transoral endoscopic insertion of foam drain

- a) Intraluminal CNP endo therapy
- b) Intracavitary CNP endo therapy
- 11) Placement of the drainage tube via the nose
- 12) Removal of the foam drain
- 13) Important notes / precautions
- 14) Reporting obligation

1. Kit description and composition

The Suprasorb® CNP endo Kit Oral consists of the following components:

- 1 **1 overtube** with grip tab (polyvinyl chloride with inner metal coil, 32 cm)
- 1 **nasal tube** (polyvinyl chloride, 50 cm) with atraumatic rounded end
- 1 **foam drain (drain)** made of an open-pore foam (polyurethane, 20 mm diameter; 100 mm length) and a transparent drainage tube (polyurethane, 2.8 mm inner diameter and 4.4 mm outer diameter, 120 cm length) with proximal retention seam and distal loop
- 1 **tube clamp** (polypolypropylene)
- 1 **nasal dressing** (nylon, polyurethane, polyethylene, medical adhesive) 83 x 70 mm

2. Required material (not included in Kit)

- Obligatory:
 - Suprasorb® CNP endo therapy unit with a connected exudate chamber
 - 1 flexible endoscope (we recommend a standard gastroscopie with a diameter of less than 10 mm)
 - Endoscope tower
 - CO₂ gas generator
 - Endoscopic grasper (non-cutting and with no barbs that could get caught in the foam)
 - Endoscopic lubricant
 - Bite block
 - Sterile scissors
- Optional:
 - Sterile needle holder and sutures
 - Narrow lumen nasal endoscope
 - Magill forceps
 - Intubation spatula

3. Intended purpose

Translumenal CNP endo therapy unit for negative pressure therapy in the oesophagus to support defect and wound healing.

4. Indications

- defects, injuries and wounds in the oesophagus
- perforations (iatrogenic perforations, spontaneous perforations)
- anastomoses/insufficiencies

5. Contraindications

- defects with involvement of the bronchial system (bronchus/trachea/lungs)
- direct contact of the foam drain with oesophageal varices
- pre-existing coagulation disorders (risk of bleeding)
- direct contact of the foam drain with large blood vessels
- irreversible ischaemia in the area being treated
- known intolerance or allergy to one or more components of the product
- injuries or constrictions in the mouth or throat that do not allow safe insertion of the CNP endo Overtube
- injuries or constrictions in the nose that do not allow safe insertion of the CNP endo Nasal Tube (only for placement via the nose)

6. Side effects

- bleeding after removal of the foam drain
- damage to the mucosa due to difficult removal of foam
- pain and/or discomfort after removal of the foam drain
- nose and throat soreness
- dysphagia
- bronchooesophageal fistula
- disturbance of the structure and stenosis
- pneumomediastinum
- haemorrhage
- fungal endophthalmitis (reported in isolated cases in the literature)

7. Interactions

There are no known interactions to date.

8. Storage

Store the Suprasorb® CNP endo Kit Oral in a dry place at a maximum of 25°C and protected from dust and sunlight.

9. Preparation

- When using the intracavitary or intraluminal negative pressure therapy is indicated before starting the therapy. Determine the location and size of the defect. The location of the defect determines the penetration depth of the foam drain. Use the markings of the endoscope as a guide for this and make a note of the penetration depth.
- Use CO₂ insufflation during the entire endoscopy until the therapy starts in accordance with the applicable standards.
- Use sufficient endoscopic lubricant during all placement manoeuvres to ensure the sliding ability of overtube, foam drain, grasper, nasal tube and endoscope.

Trimming the foam drain for intracavitary application:

- If needed, trim the foam using sterile scissors to the required dimension matched to the wound situation. Trimming should allow the foam to be placed fully in the wound cavity, however, the foam must not be shorter than 2 cm from the proximal retention seam, as otherwise, the distribution of negative pressure can no longer be guaranteed (see Fig. 1).
- Do not cut the foam in the immediate vicinity of the patient and remove loose foam particles before application. Then, trim a few millimetres off the drainage tube separately, so that the foam hangs over the end of the tube (see Fig. 2).
- By trimming, you automatically remove the thread loop at the end of the foam drain (into which you hook the grasper when inserting into the oesophagus). If the foam drain has been shortened, the thread loop must be replaced onto the foam or grasp the foam directly with the grasper, do not grasp the edge, instead grasp as close to the centre as possible to avoid tearing.

10. Transoral endoscopic insertion of foam drain

10.1 Intraluminal CNP endo therapy

Use a bite block when performing the therapy (see Fig. 3). In intraluminal CNP endo therapy, the open-pore gelled foam drain is inserted directly in the oesophagus, untrimmed, at the height of the defect, which the foam covers proximally and distally.

10a.1. Placement of the overtube

Apply lubricant to the inside and outside of the overtube. Then slide it over the endoscope. Insert the endoscope into the oesophagus and then, as a guide aid, insert the overtube on the endoscope until the non-reinforced end part (grip tab and funnel) of the overtube fully protrudes from the patient's mouth (see Fig. 4). Now remove the endoscope.

10a.2. Grasping the foam drain

Insert a suitable endoscopic grasper (non-cutting and with no barbs that could get caught in the foam) into the working channel of the endoscope. Then, grasp the thread loop at the end of the foam drain with the grasper. Place the endoscope behind the foam (see Fig. 5a). If you grasp the foam directly with the grasper, do not grasp the edge, instead grasp as close to the centre as possible to avoid tearing.

10a.3. Insertion of the foam drain into the oesophagus

Insert the endoscope with foam drain grasped through the overtube into the oesophagus (see Fig. 6a). While doing so, hold the overtube in place in the oesophagus with the grip tab. When repositioning the endoscope, the grasper must be held in place (e.g. by securing it to the handle of the endoscope with the small and ring finger of the left hand). Guide the endoscope to the front end of the foam drain. This allows the foam to be placed below the defect under endoscopic monitoring. Check the penetration depth of the foam drain using the markings on the endoscope (see Fig. 7a).

10a.4. Removal of the overtube

Hold the foam drain and endoscope in position while pulling the overtube out of the mouth and back on the endoscope (see Fig. 8a, step 1). Pull the drainage tube out of the overtube and connect it to the Suprasorb® CNP endo therapy unit (see Fig. 8a, step 2). Push the overtube to the proximal end of the endoscope and leave it there.

10a.5. Placement of the foam drain on the defect

Ensure that you have an endoscopic view of the defect and that the foam is placed below the defect.

Now open and remove the grasper. Under endoscopic visual control, pull the foam drain on the drainage tube back to the desired location at the level of the defect (see Fig. 9a).

In intraluminal CNP endo therapy, optimal positioning is achieved when the defect is positioned in the middle of the foam of the foam drain.

10a.6. Starting the therapy

Start to establish negative pressure. We recommend a negative pressure of -125 mmHg. The negative pressure can be adjusted according to the wound situation. The foam attaches to the tissue

through suction and the oesophageal lumen collapses. Now irrigate with an irrigation fluid. The therapy is working when air bubbles and fluid rise in the drainage tube.

Only remove the endoscope after checking the position of the foam drain a final time. The x-ray detectability of the drainage tube allows the subsequent use of x-ray to check the positioning.

10b) Intracavitary CNP endo therapy

Use a bite block when performing the therapy (see Fig. 3). In intracavitary endoscopic CNP endo therapy, the gelled foam drain is placed in the extraluminal wound cavity through a defect. The foam drain should be trimmed depending on the size of the wound cavity. Proceed according to Chapter 9) Preparation in the section "Trimming the foam drain for intracavitary application".

10b.1. Placement of the overtube

Apply lubricant to the inside and outside of the overtube. Then slide it over the endoscope. Insert the endoscope into the oesophagus and then, as a guide aid, insert the overtube on the endoscope into the oesophagus until the non-reinforced end part (grip tab and funnel) of the overtube fully protrudes from the patient's mouth (see Fig. 4). Now remove the endoscope.

10b.2. Grasping the foam drain

Insert a suitable endoscopic grasper (non-cutting and with no barbs that could get caught in the foam) into the working channel of the endoscope. Then, grasp the thread loop at the end of the foam drain with the grasper. If the foam is trimmed, grasp the centre of the foam. Place the endoscope behind the foam (see Fig. 5b). If you grasp the foam directly with the grasper, do not grasp the edge, instead grasp as close to the centre as possible to avoid tearing.

10b.3. Insertion of the foam drain into the oesophagus

Insert the endoscope with foam drain grasped through the overtube into the oesophagus (see Fig. 6b). While doing so, hold the overtube in place in the oesophagus with the grip tab. When repositioning the endoscope, the grasper must be held in place (e.g. by securing it to the handle of the endoscope with the small and ring finger of the left hand). Check the penetration depth of the foam drain using the markings on the endoscope (see Fig. 7b).

10b.4. Removal of the overtube

Hold the foam drain and endoscope in position while pulling the overtube out of the mouth and back on the endoscope (see Fig. 8b, step 1). Pull the drainage tube out of the overtube and connect it to the Suprasorb® CNP endo therapy unit (see Fig. 8b, step 2). Push the overtube to the proximal end of the endoscope and leave it there.

10b.5. Placement of the foam drain in the wound cavity

Guide the endoscope to the front end of the foam drain. Place the foam, under endoscopic visual control, through the defect and in the extraluminal wound cavity. Fully insert the foam into the wound cavity (see Fig. 9b). If you discover a narrow wound cavity during the intracavitary CNP endo therapy, a bougienage of the defect may be necessary if the foam drain cannot be inserted into the wound cavity due to its size. If this procedure is not possible, change to intraluminal endoscopic CNP endo therapy if appropriate. Now open and remove the grasper. Then, pull the endoscope back until you are looking at the defect opening from the oesophageal lumen.

In intracavitary CNP endo therapy, optimal positioning is achieved if the inserted foam of the foam drain is fully enveloped by the wound cavity.

10b.6. Starting the therapy

Start to establish negative pressure. We recommend a negative pressure of -125 mmHg. The negative pressure can be adjusted according to the wound situation. The foam attaches to the tissue through suction and the wound cavity begins to collapse. Now irrigate with an irrigation fluid. The therapy is working when air bubbles and fluid rise in the drainage tube and the opening of the defect visibly collapses around the drainage tube.

Only remove the endoscope after checking the position of the foam drain a final time. The x-ray detectability of the drainage tube allows the subsequent use of x-ray to check the positioning.

11. Placement of the drainage tube via the nose (optional at the discretion of the treating physician)

- 11.1. Insert the nasal tube through a nostril and slightly out of the mouth. Take care to ensure that the rounded tip of the nasal tube is inserted first (see Fig. 10)! Use Magill forceps and an intubation spatula as needed. Under no circumstances should the nasal tube be inserted with the open side first, as indicated by a yellow marking on the tube. If placement via the selected nostril causes problems, stop the procedure and perform the placement via the other nostril.
- 11.2. Then, cut the rounded tip of the nasal tube using sterile scissors to allow it to be subsequently connected to the drainage tube (see Fig. 10). Do not cut the tip of the nasal tube off over the patient's mouth in case it falls in.
- 11.3. Clamp the drainage tube with the supplied tube clamp 5–10 cm in front of the mouth so that the negative pressure does not leak.
- 11.4. Pause the Suprasorb® CNP endo therapy unit and disconnect the drainage tube (the connector remains at the exudate chamber tube).
- 11.5. Insert and connect the drainage tube with the nasal tube (see Fig. 11).
- 11.6. Pull the drainage tube out through the nose with the nasal tube. Ensure that the tube clamp is not dragged into the mouth!
- 11.7. Disconnect the nasal and drainage tube. Discard the nasal tube (see Fig. 12, step 1).
- 11.8. Reconnect the drainage tube to the Suprasorb® CNP endo therapy unit (see Fig. 12, step 2).
- 11.9. The drainage tube is now inserted into the mouth and the drainage tube is inserted into the mouth.
- 11.10. Open and remove the tube clamp.
- 11.11. Now pull the drainage tube further out of the nose until it is straightened in the throat (see Fig. 13).
- 11.12. Check the position of the foam drain a final time! Check the position (digitally or using an endoscope) of the drainage tube in the throat. The tube should be positioned in the middle of the defect.
- 11.13. Fix the drainage tube with the nasal dressing:
 - Clean, remove excess oil and dry the nasal bridge.
 - Remove the dressing's protective film and place it on the nasal bridge, with the extension of the dressing protruding over the tip of the nose. Wait a few seconds after applying the dressing so that full adhesion can develop (see Fig. 14).
 - Then, wrap the lower end of the dressing around the drainage tube in a spiral pattern towards the tip of the nose (see Fig. 14).
 - Do not use in case of known allergies to wound dressing adhesives, skin defects in the area of application or on heavily sweating patients. Do not allow the adhesive elements to come into contact with mucous membranes.
 - Mark the drainage tube at the body exit site. The marking assists you in detecting slippage of the foam drain during the therapy.

12. Removal of the foam drain

Before removal, endoscopically check that the position of the foam drain is correct.

Pause the therapy and disconnect the drainage tube (the connector remains on the exudate chamber tube) from the Suprasorb® CNP endo therapy unit.

Carefully roll the nasal dressing from the nose downwards until it is completely detached from the patient's skin. Secure the drainage tube in place and carefully roll the dressing off it. Do not use dressing removers, benzene, creams or oils.

Pull on the drainage tube carefully and continuously. A firm adhesion of the tissue to the foam due to the negative pressure exerted is a sign that the therapy is working properly; often the foam only detaches gradually from the suctioned tissue. Very firmly suctioned tissue can be detached from the foam with careful sideways movements of the blunt tip of the endoscope and with the aid of endoscopic instruments. If this does not help, secure the drainage tube with an adhesive strip and pause the therapy for up to two hours. You can then attempt to detach the foam drain again.

If the foam drain can be detached, first cut off the drainage tube near the nose and then grasp the tube in the throat and reverse placement via the nose. If necessary, use Magill forceps and an intubation spatula to help you.

Remove the foam drain under endoscopic monitoring. Position the endoscope with a view of the upper end of the foam.

Full removal of the foam drain through the patient's mouth. Do not remove the foam drain through the nose!

Using the endoscope, check and record that the foam has been fully removed. If residual foam remains in the wound, it must be removed using an endoscopic grasper.

13. Important notes / precautions

- Training is required to use this Suprasorb® CNP endo Kit. You can receive this training from qualified personnel at L&R or through an authorised L&R distribution partner.
- Before using the Kit, ensure that the components are complete and intact.
- In the implementation of CNP endo therapy, the components of the Suprasorb® CNP endo Kit may only be replaced by components of other Suprasorb® CNP endo Kits.
- The implementation of the endoscopic procedure must be assessed by the respective healthcare professional responsible for the therapy, taking into consideration the benefits/risks for the patient.
- The CNP endo therapy may only be performed by physicians who are experienced in interventional flexible endoscopy of the gastrointestinal tract.
- The overtube, endoscope and foam drain must be inserted and removed carefully and never suddenly or with excessive force.
- If a severe bleeding complication is suspected (e.g. in case of abnormal drainage secretions), immediately perform endoscopic monitoring of the therapy under surgical conditions and in readiness to perform surgery.
- The user decides when the therapy is to be ended depending on the progress of wound healing. We recommend replacing the foam after therapy intervals of 48 to no more than 72 hours.
- Perform early endoscopic monitoring of the therapy in case of suspected therapy failure, for example, due to displacement or obstruction of the foam.
- The overtube contains metallic elements and must not be used around an MRI (see Chapter 1).

14. Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

15. Reporting obligation

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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