Dauerbinde® Perfekta®

Langzug-Binde

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Long stretch bandage

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Ligadura de alta elasticidade

Benda ad allungamento elevato

Lange-rekzwachtel

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Elastinen kompressioside

Obinadlo s dlouhým tahom

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Dugoelastični zavoj

Pika venivusega side

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en Instructions for Use

Product description and performance characteristics
Dauerbinde® K / Perfekta® kräftig is a long stretch bandage with approx. 175% extensibility.
Dauerbinde® F / Perfekta® super is a long stretch bandage with approx. 160% extensibility.

Product composition
Dauerbinde® K / Perfekta® kräftig: 58% cotton, 24% polyamide, 8% elastane, 10% viscose
Dauerbinde® F / Perfekta® super: 63% cotton, 28% polyamide, 9% elastane

Intended purpose
For compression, fixation of dressings, support and relief and for immobilisation of limbs.

Indications
phlebological and lymphological indications: varicosis, deep venous insufficiency, thrombophlebitis, phlebothrombosis, post-thrombotic syndrome, stage I to III b chronic venous insufficiency according to Widmer, venous and mixed ulcers, oedema (phleboedema, lymphoedema, lipoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery.
contusions, sprains, tendon injuries and joint dislocations

Contraindications

Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects

In very rare cases skin irritations and / or allergies may occur.

Warnings and precautions

- Due to high resting pressure, long stretch bandages must be removed at night.
- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Slight bluish tinge of the toes is normal, this should disappear as soon as the patient starts to walk.
- If the symptoms do not disappear (pain and bluish tinge), the bandage must be removed and applied again with less compression.
- The patient must inform their doctor if sharply increasing pain occurs while

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