Gazin® RK

Spitztupfer mit Röntgenkontrast, unsteril
Tampons pointus avec contraste à la radiographie, non stériles
Pointed swabs RK, non-sterile
Tampones de gasa con medio de contraste, no estéril
Tamponi con punta affusolata con contrasto radiopa-co, non sterili
Puntdeppers met röntgencontrast, niet-steriel
Hegyes törlő, röntgenkontrasztot adó, nem steril

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en Instructions for Use

Product description
Gazin® pointed swabs RK, non-sterile, consist of 20 thread count cotton gauze compliant with EN 14079. The X-ray detectable thread contained in the swabs consists of barium sulphate filaments and is wound with polypropylene.

Product composition
100% cotton, X-ray detectable thread (blue, polypropylene, polyester thread, barium sulphate).

Intended purpose
Gazin® pointed swabs are used for absorbing body fluids, for haemostasis or as tamponade in the course of surgical procedures.

Usage information/Precautions
This product is non-sterile. Gazin® pointed swabs must be sterilised prior to surgically invasive use; resterilisation is not permitted. The product is intended for single use only.

Sterilisation instructions
1. Device requirements
Only sterilizers, which comply to the requirements of the norms EN 285 or EN 13060 and which use a pre-vacuum system for the sterilisation are allowed to be used. The use of the gravitational method is not allowed in the EU.

2. Sterilisation preparation
Before the product is sterilised, a suitable sterile barrier system must be selected, validated and used by the person responsible for the sterilisation (sterilising technician) in accordance with ISO 11607. The sterile barrier must be big enough for the product so that the seal is not under tension. The retail packaging is unsuitable for sterilisation.

3. Sterilisation
The sterilisation procedure for the product must be validated by means of moist heat at 134°C for 7 minutes in accordance with ISO 17665. Do not use after sterilisation if package is damaged or accidentally opened.

4. Storage
After removing the required quantity of Gazin® pointed swabs, non-sterile, seal the packaging and store dry and protected from dust. Sterilised products must be stored dry and protected from dust in microbe-resistant sterile packaging.

Additional information
The above instructions were assessed by Lohmann & Rauscher as suitable for preparation for sterilisation. The sterilising technician is responsible for ensuring that the sterilisation carried out with the employed equipment, materials and staff achieves the required results. To this end, validation of the sterilisation process and routine monitoring of the procedure is normally required. Deviations from the procedure are the responsibility of the sterilising technician. National requirements for steam sterilization must be considered by the user.

Disposal
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

General instructions
In case of a serious incident, contact the manufacturer and the responsible health authority.

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