Gebrauchsanweisung
Mode d'emploi
Instructions for use
Instrucciones de uso
Istruzioni per l'uso
Gebruiksaanwijzing
Návod k použití
Navodila za uporabo

Suprasorb® CNP P3

Suprasorb®

Controlled
Negative
Pressure

Inhalt

Deutsch	5 - 56
Français	57 - 108
English	109 - 160
Español	161 - 212
Italiano	213 - 264
Nederlands	265 - 316
Čeština	317 - 368
Slovenščina	369 - 420

Suprasorb® CNP P3

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Suprasorb®

Controlled
Negative
Pressure

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Suprasorb®

Controlled
Negative
Pressure

The safety measures of the **Suprasorb® CNP P3** therapy unit comply with the regulations of the **German Medical Devices Act (MPG)** and the generally recognised rules of technology.

Suprasorb® CNP P3 is approved in accordance with EU Council Directive 93/42/EEC on medical devices and fulfils in particular the essential requirements of Annex I of the same directive.

Suprasorb® CNP P3 is a medical suction device and has been classified as a Class IIa medical device in accordance with Annex IX, Rule 11 of EC Directive 93/42/EEC.

Suprasorb® CNP P3 satisfies the requirements of the IEC 60601-1-2/EN 60601-1-2 “Electromagnetic Compatibility – Medical electrical equipment” standard.

This minimises electromagnetic interference and interactions.

The **quality management system** in place at L&R has been certified in accordance with the pertinent international standards.

Last revised 07-2017

Copyright © 2017

© This document, including excerpts, must not be republished.
The contents as a whole or parts thereof may not be reproduced,
copied or provided to third parties in any shape or form without
the written consent of L&R.

Subject to change without notice. Errors and omissions excepted.

TABLE OF CONTENTS

1	Information for the user	114
1.1	Use of the instructions for use	114
1.2	Symbols and their meaning	115
1.3	Glossary of definitions	117
1.4	Intended purpose	118
1.4.1	Indications	118
1.4.2	Contraindications	118
1.4.3	Precautionary measures	119
1.4.4	Restrictions on use	119
1.5	Information for the user	119
1.6	If the patient is discharged	120
2	Information on Suprasorb® CNP P3 and its components	120
2.1	Suprasorb® CNP P3 scope of delivery	121
2.2	Suprasorb® CNP P3 properties	122
2.3	Container system	122
2.4	Filter system	123
2.5	Gelling agent	123
2.6	Battery	123
2.7	Information on charging the battery	124
2.8	Pressure settings	124
2.9	Warranty	125
3	Preparation for commissioning	126
3.1	Prior to commissioning	126
3.2	Function test	127
3.3	Insertion of the Suprasorb® CNP Exudate Pouch and connection of the container system	128
3.4	Changing the Suprasorb® CNP Exudate Pouch	130
3.5	Positioning the Suprasorb® CNP P3	132
3.6	Fixing the Suprasorb® CNP P3 in the bed mount/universal mount	134
3.7	Removing the Suprasorb® CNP P3 from the bed mount/universal mount	135
3.8	Connection to wound drainage system and power supply	135

4	Operation by qualified healthcare professionals	137
4.1	Switching the device on and off	137
4.2	Screen saver and screen lock	137
4.3	Start screen	138
4.4	Preset therapy mode	138
4.5	Changing settings	138
4.6	Setting the therapy mode	138
4.7	Setting the sensitivity	139
4.8	Confirming the therapy setting	140
4.9	Starting the therapy and seal level	141
4.10	Adjusting the language, volume, brightness, date and time	142
4.11	Viewing the chronology of the changed settings	143
5	Operation by patient	144
5.1	Switching the device on and off	144
5.2	Screen saver and screen lock	144
5.3	Start screen	145
5.4	Preset therapy mode	145
5.5	Starting the therapy and switching to Pause mode	145
5.6	Adjusting the language, volume, brightness, date and time	146
5.7	If the doctor has granted permission to shower	147
6	Handling error messages	148
6.1	Function test	148
6.2	Device does not start	148
6.3	Error messages on the screen and solutions	148
7	Cleaning of Suprasorb® CNP P3 therapy unit	151
7.1	Cleaning and disinfection of the device surfaces and the container	151
7.2	Maintenance and service	152
8	Transport, storage and disposal	153
8.1	Decontamination before shipping	153
8.2	Storage	153
8.3	Disposal	154
9	Technical data of the Suprasorb® CNP P3 therapy unit and EMC information	155
10	Ordering options	160
11	Contact information and hotlines	160

1 Information for the user

1.1 Use of the instructions for use

Read the instructions for use carefully before using the **Suprasorb® CNP P3** for the first time.

Read Section 3 “Preparation for commissioning” before commissioning the device.

Sections 4 and 5 are intended for different user groups:

- **If you are a qualified healthcare professional** Read Section 4 “Operation of the **Suprasorb® CNP P3** by qualified healthcare professionals” for information on operation of the device.
- **If you are the patient or a friend/relative** Read Section 5 “Operation of the **Suprasorb® CNP P3** by the patient” for information on operation of the device.

Pay particular attention to the safety information at the beginning of each section to avoid risks and hazards.

Keep these instructions for use in a safe and readily accessible place.

If passing on the **Suprasorb® CNP P3** therapy unit to third parties, always include the instructions for use.

Quick guide to these instructions for use:

- To insert or replace the exudate pouch:
Go to Section 3.3. and 3.4.
- To position the device safely:
Go to Section 3.5.
- To use the bed mount or universal mount:
Go to Section 3.6.
- To connect the device to the power supply:
Go to Section 3.8.
- To handle any error messages displayed on the screen:
Go to Section 6.3.

1.2 Symbols and their meaning

General:

Symbol	Meaning	Symbol	Meaning
 WARNING	Risk of serious physical injury	 CAUTION	Risk of physical injury and harm to health
 PLEASE NOTE	Other useful information/ advice on how to avoid material damage		

On the device and packaging:

Symbol	Meaning	Symbol	Meaning
	Avoid wetting		Do not dispose of this device in household waste.
	Protection class II		Order number
	Atmospheric pressure range		Serial number
	Humidity range		Lot number
	Follow the instructions for use		Date of manufacture
	Type BF applied part		Manufacturer
	Temperature range		Do not use if packaging not intact!
	Power supply		Not for reuse
IP23	International protection / degree of protection The degree of protection gives the device's protection against contact and ingress of liquids. The Suprasorb® CNP P3 therapy unit is protected against the insertion of fingers and splashing water at an inclination of up to 60°.		

On the screen:

Symbol	Meaning
	Key for switching the Suprasorb® CNP P3 on and off
	Power supply is connected when green dot is illuminated
	Browse forwards/backwards through the list of possible settings
	Continuous therapy mode
	Intermittent therapy mode
	Pressure
	Time
	Decrease/increase value
	OK (Enter)
	Cancel
	Battery full
	Battery low
	Key lock: Activates automatically during operation and can be unlocked by holding the key down
	General settings
	Back
	Start therapy
	Pause therapy
	Setting the therapy mode
	Therapy sensitivity
	Language selection
	Signal volume
	Log file
	Setting the time and date
	Setting the screen brightness
	Mute audible alarms for a period of two minutes

1.3 Glossary of definitions

C

Collapse	Cave in on oneself / crumple. Collapsing of the exudate pouch is undesirable. In such cases, it is important to establish whether it has been inserted correctly in the device or if any leaks can be found. In contrast, collapsing of the wound dressing indicates that negative pressure has been established in the wound and that the therapy is functioning as intended.
Contamination	Contamination is when micro-organisms such as bacteria and fungi as well as viruses colonise a surface, compromising its cleanliness.

E

EMC	Electromagnetic compatibility: The device's capacity not to interfere with other devices as a result of adverse electrical and/or electromagnetic effects or to be affected by other devices.
Exudate	The term exudate refers to all accumulating fluids and particles which can develop or be present in a wound. The exudate is suctioned out of the wound with the help of the Suprasorb® CNP P3 therapy unit and collected in an exudate pouch.

H

Haemostasis	Haemostasis is the sum of the physical processes which combine to stop bleeding.
-------------	--

L

L&R	L&R is the abbreviation used in the instructions for use for Lohmann & Rauscher International GmbH & Co. KG with its headquarters in Rengsdorf, Germany.
-----	--

O

Oversuction	Oversuction is when exudate is sucked into the interior of the device.
-------------	--

R

Reprocessing	Reprocessing is required whenever there is a change of patient. Reprocessing means that all the single-use products are disposed of and the Suprasorb® CNP P3 therapy unit wipe-down disinfected. The Suprasorb® CNP P3 therapy unit must only be opened by L&R or an authorised L&R service partner.
--------------	--

S

Sensitivity	The device's sensitivity specifies how sensitively the device is set and how quickly alarms are triggered. The higher the sensitivity setting, the more sensitively the device reacts and the quicker the alarm is triggered.
-------------	---

T

Therapy unit	Designates the Suprasorb® CNP P3 medical suction device incl. power supply and two containers – Suprasorb® CNP Container 250 ml and Suprasorb® CNP Container 1,000 ml. They are employed during Suprasorb® CNP (Controlled Negative Pressure) therapy.
--------------	---

1.4 Intended purpose

The **Suprasorb® CNP P3** therapy unit is designed to drain exudate by suction during Suprasorb® CNP (Controlled Negative Pressure) therapy. It is intended for patients likely to benefit from negative pressure wound care.

1.4.1 Indications

The **Suprasorb® CNP P3** therapy unit can be used on both uninfected and contaminated, colonised, critically colonised or infected wounds with light, moderate or heavy amounts of exudate and which are superficial or deep.

Possible areas of application:

- lower leg ulcers (venous, arterial, mixed)
- decubitus
- diabetic foot lesions
- post-traumatic and post-operative wounds
- initially infected wounds following surgical debridement
- soft tissue injuries (lacerated contusions, degloving)
- injuries exposing the underlying bone or bradytrophic tissue
- wounds following splitting of a compartment on an extremity
- conditions following skin grafts (e.g., mesh graft, up to reliable integration of the graft on day 5/6)
- sternal wound infections following surgical debridement
- open abdominal procedures including treatment of fistulae (in combination with Suprasorb® CNP Drainage Film only)
- treatment of 1st and 2nd degree burns

1.4.2 Contraindications

Suprasorb® CNP P3 is contraindicated for the following applications:

- exposed vessels likely to be compromised by negative pressure; the same applies for vascular anastomoses
- coagulation disorders (risk of haemorrhaging)
- slight to heavy acute bleeding in the wound area following injury/surgical debridement
- necrotic tissue
- unexplored fistulas
- untreated osteomyelitis
- malignant wounds
- exposed organs without suitable organ protection
- dry wound conditions
- 3rd degree burns
- use on the sternum
(with the exception of sternal wound infections following surgical debridement)

1.4.3 Precautionary measures

Precautionary measures need to be taken for:

- patients under anticoagulation therapy and/or with active bleeding
- patients with difficult wound haemostasis
- application of the system in the immediate vicinity of blood vessels, organs, muscles and fascias and bradytrophic tissue (tendons, ligaments or cartilage)
- irradiated vessels and tissues
- bone fragments with sharp edges
- uncooperative patients
- wound cavities and pockets with hidden areas
- patients with ischaemia and/or circulatory disorders; the attending doctor must check whether Suprasorb® CNP therapy is applicable or whether it may be necessary to implement other therapy measures (e.g., revascularisation) first

1.4.4 Restrictions on use

- medical facilities where potential equalisation is required (e.g., cardiac surgery)
- potentially explosive atmospheres
- aircraft
- areas with electromagnetic radiation, e.g., MRI suites

1.5 Information for the user

- Use the **Suprasorb® CNP P3** in accordance with its intended purpose.
- Use the **Suprasorb® CNP P3** in accordance with its instructions for use.
- The **Suprasorb® CNP P3** and its components should always undergo a visual inspection and a functional check prior to use.
- Do not allow unauthorised persons to open the **Suprasorb® CNP P3** therapy unit.
- Assembly, setup, add-ons, maintenance and repairs must be performed by authorised personnel.
- The **Suprasorb® CNP P3** must undergo a check in accordance with DIN EN 62353 performed by L&R or an L&R authorised service partner every 18 months (see Section 7.2).
- The healthcare professionals are trained in the use of the device by L&R or an authorised L&R sales partner. The training sessions are performed once with reference to the instructions for use – Sections 2, 3, 4, 6 and 11 for the qualified healthcare professionals and Sections 2, 3, 5, 6 and 11 for the patient – and documented in the medical devices book.
- The **Suprasorb® CNP P3** has a product service life of 3 years.

1.6 If the patient is discharged


WARNING

Important points to follow when performing mobile CNP at home
Patients need to be instructed in the proper use of **Suprasorb® CNP P3** if therapy is to be performed in their home environment. In addition, it is imperative that trained healthcare professionals monitor the wound dressing and the negative pressure system regularly.

- Use of the Suprasorb CNP therapy is suitable only for patients whose clinical condition permits treatment at home (adequate haemostasis, low risk of haemorrhaging).
- The patient, relatives or carer must be capable of reading and understanding the safety information, reacting appropriately to error messages and following the instructions.
- **The Suprasorb® CNP Container 1,000 ml is not intended for treatment at home.**

2 Information on Suprasorb® CNP P3 and its components


Suprasorb® CNP P3

- A** Suprasorb® CNP Container 250 ml
(the Suprasorb® CNP Exudate Pouch 250 ml is inserted in this)
- B** Exudate pouch tube plus tube clamp as part of the exudate pouch
- C** Screen
- D** Control panel below screen, alongside power indicator


 Power on and off key

-  Mute key for audible alarms
-  Control LED for connected power supply
- E Power supply port
- F LED bar:

The LED bar is green when the therapy is running without any problems.



The LED bar goes orange when an error message requires attention.

2.1 Suprasorb® CNP P3 scope of delivery

The case contains:

- A 1 Suprasorb® CNP P3 medical suction device
- B 1 Suprasorb® CNP Container 250 ml
- C 1 Suprasorb® CNP Container 1,000 ml
- D 1 power supply TR15RA150-E-160E24 incl. European adapter
- E 1 Suprasorb® CNP Bed Mount
- F 1 Suprasorb® CNP Carry Bag
- G 1 instructions for use for qualified healthcare professionals and patient, charging information, medical devices book

2.2 Suprasorb® CNP P3 properties

CAUTION

Use of the carry bag

- The Suprasorb® CNP Carry Bag for mobile use supplied with the **Suprasorb® CNP P3** is intended for single use per patient. It should not be used for other patients but be disposed of once used.

The **Suprasorb® CNP P3** is a pump designed for stationary and mobile use for the medical drainage of exudate during Suprasorb® CNP therapy. It is used for drainage in the low negative pressure range and can be employed in both hospitals and practices, as well as when transporting patients and in home care.

The **Suprasorb® CNP P3** is rechargeable and easily portable. The **Suprasorb® CNP P3** therapy unit operates on the internal rechargeable battery or the supplied power supply, which can also be used to charge the battery.

The negative pressure is generated by a maintenance-free membrane pump with an electromechanical motor. Once turned on, the negative pressure pump generates negative pressure in the tubing system and container to suction the exudate.

The exudate is routed away from the patient and collected in the container system. The integrated overflow prevention system triggers the “System closed” alert and stops the therapy as soon as the container system is full.

The **Suprasorb® CNP P3** must only be operated with the supplied container system.

The expected product service life of the **Suprasorb® CNP P3** is 3 years.

2.3 Container system

The **Suprasorb® CNP P3**'s container system is composed of a container (250 ml or 1,000 ml) and an exudate pouch, which is inserted in the container in the corresponding size (250 ml or 1,000 ml). The user has the choice of which container system to use depending on the wound to be treated and the amount of exudate.

PLEASE NOTE

The Suprasorb® CNP Container 250 ml is used in combination with the Suprasorb® CNP Exudate Pouch 250 ml. The Suprasorb® CNP Container 1,000 ml is used in combination with the Suprasorb® CNP Exudate Pouch 1,000 ml.

2.4 Filter system

There are a total of three filters integrated in the **Suprasorb® CNP P3**: An internal filter inside the **Suprasorb® CNP P3** device and two external filters in the exudate pouch.

The external charcoal/bacterial filter in the exudate pouch prevents oversuction. Draining stops when the exudate pouch needs to be changed. The second filter is a ventilation filter.

The filter inside the **Suprasorb® CNP P3** is also a charcoal/bacterial filter and protects the device against oversuction and contamination.

Should the internal filter become blocked, the **Suprasorb® CNP P3** must be properly reprocessed by the L&R service centre or an L&R authorised service partner.

The activated charcoal filter in the container prevents the development of odours.

2.5 Gelling agent

The exudate drained into the exudate pouch thickens within a few minutes due to the gelling agent contained in the exudate pouch. To avoid leaks – for example, when changing the exudate pouch – the tube clamp on the exudate pouch tube should be closed off.

The exudate pouch incl. exudate pouch tube and tube clamp is intended for single use only.

PLEASE NOTE

The internal charcoal/bacterial filter is not intended for reuse and must be replaced if it comes into contact with the exudate (blockage), once the filter service life has elapsed or during servicing/repair work performed by L&R or an L&R authorised service partner.

2.6 Battery

Battery, rechargeable	10.8 V, 2.7 Ah – lithium ion
Charging time when battery is empty	6-7 hours
Charging time with approx. 50% residual charge	3-3.5 hours
Operating time	With battery: 10-20 hours

2.7 Information on charging the battery

On delivery and prior to use of the **Suprasorb® CNP P3** therapy unit, it is expressly recommended to fully charge the battery and recharge after the first few uses. This allows achievement of the maximum number of approximately 300 charging cycles.

Suprasorb® CNP P3 is equipped with a lithium ion battery, which, unlike standard rechargeable batteries, only exhibits a very low, temperature-dependent self-discharge of 3-10% per month. The **Suprasorb® CNP P3** should be stored and charged taking the ambient conditions specified in the technical data into account, ideally at room temperature.

The **Suprasorb® CNP P3** including battery must **never** be stored uncharged!

Lithium ion batteries do not have a memory effect. They can and thus should be **recharged at any time** following a successful initial charge (see above).

However, **avoid charging the battery frequently for short periods of time.**

Even though the **Suprasorb® CNP P3** battery includes safety precautions to protect against deep discharging, the charging information above should still be taken into consideration. The battery is also protected against overheating during charging. Should the battery temperature be exceeded during charging due to ambient conditions other than intended, charging will be temporarily interrupted to allow cooling. This is a safety precaution to protect the battery.

PLEASE NOTE

The **Suprasorb® CNP P3** can be used when charging via the power supply!

2.8 Pressure settings

The pressure settings range between -20 and -200 mmHg (in 5 mmHg intervals).

Our recommendations:

- (a) Negative pressure setting in Continuous mode:
-60 to -80 mmHg
- (b) Negative pressure setting in Intermittent mode:
Max: -60 to -80 mmHg at a 2-5 minute interval
Min: -40 to -60 mmHg at a 2-5 minute interval

The negative pressure settings can be customised by the qualified healthcare professional as soon as the **Suprasorb® CNP P3** is switched on. The factory default pressure is -80 mmHg.

2.9 Warranty

The warranty for the **Suprasorb® CNP P3** runs for 2 years. This period is neither extended nor renewed when warranty work is performed. The battery comes with a 6-month warranty. The warranty does not extend to wearing parts or consumables.

L&R is only responsible for consequences for safety, reliability and specified performance capacity if:

- original accessories and spare parts are used;
- maintenance and repairs are performed by L&R or a specialist authorised by L&R;
- the device in question is operated in accordance with the instructions for use and not contrary to the intended purpose.

L&R does not guarantee the error-free functioning of the **Suprasorb® CNP P3** nor does it accept liability for any material damage or personal injuries sustained if:

- accessories and spare parts are used which are not original;
- the information in these instructions for use is not heeded;
- assembly, setup, modifications, updates or repairs are not performed by L&R or a specialist authorised by L&R;
- the security seal is removed or broken;
- another charger is used.

PLEASE NOTE

All warranty claims become void if the **Suprasorb® CNP P3** is opened by unauthorised persons, the security seal is removed/broken, repairs are performed by unauthorised persons or the device is not used in accordance with the guidelines.

3 Preparation for commissioning

3.1 Prior to commissioning

Always keep an additional exudate pouch for the **Suprasorb® CNP P3** available, as this is obligatory for safe operation!

- Read the instructions for use through completely before you use the **Suprasorb® CNP P3** therapy unit for the first time.
- Remove the **Suprasorb® CNP P3** and its accessories from the Suprasorb® CNP Case.
- Always position the **Suprasorb® CNP P3** therapy unit on a solid, flat surface, paying attention to the correct positioning of the device (horizontally or in the bed mount).
- All the tubes and the power supply need to be checked for damage each time before you use the **Suprasorb® CNP P3** therapy unit. Avoiding kinking the tubes when connecting the tubes. Ensure that the container system and tubes are connected correctly before switching on the device.
- Fully charge the battery prior to first use.
- Perform a function test (see Section 3.2).

 WARNING	Strangulation <ul style="list-style-type: none">• Ensure that the therapy is administered in such a way that the tubes and cables do not pose a strangulation risk.
 CAUTION	Setting the volume <ul style="list-style-type: none">• Before using the Suprasorb® CNP P3, check that the current volume of the alarms is appropriate for the environment. Ensure that you can hear the acoustic alarm and thus handle the error message promptly.
 PLEASE NOTE	Electromagnetic interactions <ul style="list-style-type: none">• Medical electrical devices are subject to special safety precautions in terms of electromagnetic compatibility and need to be installed and commissioned in accordance with the EMC information (see Section 9)!
 PLEASE NOTE	Incorrect accessories and spare parts <ul style="list-style-type: none">• The use of accessories and spare parts not recommended by L&R can have a detrimental effect on the safety and functioning of the device.• Damage caused by the use of accessories and spare parts other than recommended or improper use of the device is not covered by any warranty.• Use original accessories and spare parts only.
 PLEASE NOTE	Suctioning of other liquids <ul style="list-style-type: none">• Never suction flammable, corrosive or explosive liquids or gases.

 <p>PLEASE NOTE</p>	<p>If you drop the device</p> <ul style="list-style-type: none"> • Do not operate the device if the housing displays visible damage. Contact the hotline or send the device directly to the specified address.
 <p>PLEASE NOTE</p>	<p>How to avoid damage caused by liquids</p> <ul style="list-style-type: none"> • Do not use the device in rooms with high humidity, when bathing or showering. • Avoid wetting the power supply, plug and screen. • Never immerse the device in water or any other liquids (even when switched off).
 <p>PLEASE NOTE</p>	<p>Heat accumulation and dust</p> <ul style="list-style-type: none"> • The power supply and therapy unit should not be covered during operation, with the exception of storage in the carry bag provided for this purpose. • Protect the therapy unit and charging cable from direct sunlight and other sources of heat such as radiators. • Protect the device and accessories from excessive accumulation of dust.

On delivery and prior to first use of the **Suprasorb® CNP P3** therapy unit, it is expressly recommended to fully charge the battery and recharge after the first few uses.

3.2 Function test

Perform a function test on the charged therapy unit with the exudate pouch and container connected before you use the **Suprasorb® CNP P3** for therapy.



Proceed as follows:

1. Switch the **Suprasorb® CNP P3** on as described in the Section “Switching the device on and off”.
2. Start the therapy and let the therapy unit run without a connected wound dressing. The message “System is open” should appear after 30-40 seconds. If “System open” appears, proceed to 3. However, should the message “System is blocked” appear, the internal filter of the **Suprasorb® CNP P3** is blocked and needs to be replaced by an authorised L&R service partner.
3. Confirm the “System is open” message with .
4. Close the tube clamp on the exudate pouch tube. The “System is blocked” message should appear after two minutes at the latest and be confirmed with . If the message fails to appear even after multiple attempts, have the **Suprasorb® CNP P3** checked by L&R or an authorised service partner.



3.3 Insertion of the Suprasorb® CNP Exudate Pouch and connection of the container system

 <p>CAUTION</p>	<p>Exudate pouch</p> <ul style="list-style-type: none">• Ensure that the exudate pouch for the Suprasorb® CNP P3 is replaced regularly. When the exudate pouch is full, the integrated overflow prevention system is triggered. This discontinues the therapy.• Switch the device to Pause mode when replacing the exudate pouch.
 <p>PLEASE NOTE</p>	<p>Collapsing exudate pouch</p> <ul style="list-style-type: none">• A leak in the container or the cover of the exudate pouch may result in air entering the container. This, in turn, may cause the exudate pouch to collapse.• When assembling the container system, ensure that the side plate of the exudate pouch is securely connected to the container.• Ensure that all connections are secure and established correctly.• Ensure that neither the container nor the exudate pouch is damaged.• The separate instructions for use for the exudate pouch must also be followed!

The container system of the **Suprasorb® CNP P3** comprises the Suprasorb® CNP Container, the Suprasorb® CNP Exudate Pouch with the connected exudate pouch tube and a tube clamp. There are two possible combinations here:


A) Suprasorb® CNP
Exudate Pouch 250 ml
and Suprasorb® CNP
Container 250 ml


B) Suprasorb® CNP
Exudate Pouch 1,000 ml
and Suprasorb® CNP
Container 1,000 ml

Proceed as follows:

1. Remove the exudate pouch from the packaging and extend completely.


2. Insert the exudate pouch loosely into the reusable container. During insertion, ensure that the pouch film is inserted into the container completely and is not trapped at the edges.

3. Now promptly affix the container system (Suprasorb® CNP Container with inserted Suprasorb® CNP Exudate Pouch) to the main unit and as shown in the diagram until it clicks into place completely and is securely connected.

4. The connection to the wound drainage system can now be established on the exudate pouch tube. This can be done using either the Suprasorb® CNP Drainage Pad or the Suprasorb® CNP Deep Drain, each with the corresponding connector.

See also instructions for use for Suprasorb® CNP Drainage Pad; points f, g and h of the application instructions.

The three steps for connecting the tube using the single-piece connector are:

1. Insert the tube into the connector for connection.

2. The tube is inserted, but the connector is still UNLOCKED.

3. Connector is locked at both ends.

1. Insert the exudate pouch tube TWO CENTIMETRES into the connector for connection. Push a little more if necessary.
2. The exudate pouch tube is now inserted. However, the connector is still UNLOCKED. Lock the connector by pulling the sliding ends of the connector apart.
3. Figure 3 shows the connector when LOCKED.

3.4 Changing the Suprasorb® CNP Exudate Pouch

CAUTION

Avoiding infections

Infectious and pathogenic microbes in the exudate can pose health risks.

- Wear appropriate disposable gloves when **replacing the container**.
- Replace the exudate pouch including the exudate pouch tube in accordance with the respectively applicable hygiene regulations, at the latest every week, before a change of patient or when full.
- **Reprocessing is required whenever there is a change of patient.** Reprocessing means that all the single-use products are disposed of and the **Suprasorb® CNP P3** therapy unit wipe-down disinfected. The **Suprasorb® CNP P3** therapy unit must only be opened by L&R or an authorised L&R service partner.
- The exudate and all components contaminated with it must be disposed of **in accordance with the facility's policies**.


CAUTION

Exudate pouch

- Ensure that the exudate pouch for the **Suprasorb® CNP P3** is replaced regularly. When the exudate pouch is full, the integrated overflow prevention system is triggered. This discontinues the therapy.
- Switch the device to Pause mode when replacing the exudate pouch.

As soon as the exudate pouch reaches maximum capacity, the **Suprasorb® CNP P3** displays the message “Exudate pouch is full”. Please then proceed as follows:


1. Confirm the message by pressing the green tick.


2. Close the tube clamps on the exudate pouch tube and the wound drainage system to ensure that the negative pressure in the wound is maintained.


3. Switch the **Suprasorb® CNP P3** off  (Pause mode -> see Section 5.4).

4. Detach the exudate pouch tube from the wound drainage system connector.

5. Detach the container system from the remaining unit by first pushing down and then to the side, as shown in the diagram.

6. Remove the exudate pouch from the reusable container. There are two handholds integrated in the exudate pouch (250 ml or 1,000 ml) to facilitate removal (see image). These allow you to pull the exudate pouch a little out of the container before grasping by hand and removing from the container completely.

7. Dispose of the exudate pouch including the exudate pouch tube and tube clamp in accordance with the applicable regulations for household waste.
8. **Then proceed as described in Section 3.3 and insert a new exudate pouch in the reusable container.**
9. Connect the exudate pouch tube of the new exudate pouch up to the wound drainage system (see end of Section 3.3) and fix the container back in place on the unit.
10. Check whether the tube clamps are open so as to continue the therapy. Then restart the **Suprasorb® CNP P3** and continue the desired therapy.

3.5 Positioning the Suprasorb® CNP P3


PLEASE NOTE


The device should be placed horizontally.
 The device must only be positioned vertically in the bed mount or universal mount and in the carry bag.
 The device should not be positioned above the wound, but rather at the same height or below it.

The safest way to store the **Suprasorb® CNP P3** is in a horizontal position. In addition, it can be attached to the sides of the bed using the bed mount (see Section 3.6) or the universal mount.

Horizontal storage:

Storage in the bed mount/universal mount:

Storage in the carry bag:

The Suprasorb® CNP Carry Bag should be used for mobile use. However, the doctor is responsible for deciding whether the patient's condition permits mobile use.

CAUTION

Use of the carry bag

- The Suprasorb® CNP Carry Bag for mobile use supplied with the **Suprasorb® CNP P3** is intended for single use per patient. It should not be used for other patients but be disposed of once used.


PLEASE NOTE


Use of the Suprasorb® CNP P3 with the Suprasorb® CNP Carry Bag
When using the **Suprasorb® CNP P3** with the corresponding carry bag, only use the Suprasorb® CNP Container 250 ml and the Suprasorb® CNP Exudate Pouch 250 ml. Position the **Suprasorb® CNP P3** therapy unit in the Suprasorb® Carry Bag such that the screen of the device is visible through the inspection window in the carry bag.

3.6 Fixing the Suprasorb® CNP P3 in the bed mount/universal mount

The device can be safely used both in the horizontal position and when mounted in the supplied bed mount. There is also a universal mount available for the device. The procedure for inserting the device in the bed mount and the universal mount is the same.

1. Push the **Suprasorb® CNP P3** into the bed mount/universal mount from right to left. When doing so, ensure that the side with the screen is facing you.

2. The therapy unit is inserted when the feet of the device click into place in the recesses on the bottom strut of the bed mount (see arrow at bottom of diagram above for correct positioning). For final fixation, move the green lever on the top of the bed mount/universal mount to the left.

3. The therapy unit can also be inserted in the bed mount/universal mount or removed from the same when the bed mount/universal mount is attached to the sides of the patient's bed.


3.7 Removing the Suprasorb® CNP P3 from the bed mount/universal mount

1. To remove the **Suprasorb® CNP P3** from the bed mount/universal mount, move the green lever to the right and push the therapy unit gently to the right out of the bed mount.

2. The therapy unit can also be removed from the bed mount/universal mount when the latter is attached to the sides of the patient's bed or anchored on a rail.

3.8 Connection to wound drainage system and power supply

The connection to the wound drainage system can now be affixed to the exudate pouch tube. This can be done via either the Suprasorb® CNP Drainage Pad or the Suprasorb® CNP Deep Drain. The connection is established via the single-piece connector on their tube ends.

Three steps for connecting the tube using the single-piece connector: See Section 3.3.

PLEASE NOTE

The application of the wound dressing and the connection to the therapy unit must only be performed by qualified healthcare professionals. Further information on the Suprasorb® CNP Drainage Pad, Suprasorb® CNP Deep Drain and Suprasorb® CNP Y Connector can be found in the instructions for use for the respective product.

Using the power socket on the **Suprasorb® CNP P3**, connect the unit to the mains power supply for charging or operation via the power supply.

PLEASE NOTE

Incorrect voltage power supply

Improper operation causes overvoltages to develop in the device which can then pass to the user.


- Ensure that the power supply is set to a voltage of 100-240 V alternating current prior to use.
- Only operate the device with the power supply TR15RAM150-E-160E24 input 15 V DC 1 A.

If required, there is an adapter supplied in the case with which the device can be connected up to the mains power supply in your country.

Use the supplied power supply only.

Start by connecting the power supply to the power socket on the **Suprasorb® CNP P3**. Then connect the power supply to the mains power supply.

*Connection of
Suprasorb® CNP P3
therapy unit to patients
and accessories*

4 Operation by qualified healthcare professionals

4.1 Switching the device on and off

To switch the device on, press  . This key can be found below the screen on the **Suprasorb® CNP P3** (see image). To switch the device off, hold the same key down until the screen goes black.

Power supply

On/Off

Signal off

The LED of the on/off key is green when the **Suprasorb® CNP P3** is switched on. The control LED for the power supply is green when the **Suprasorb® CNP P3** is connected to the mains power supply.

4.2 Screen saver and screen lock

If the screen goes dark during therapy, it can be made visible again by simply touching it. However, the editing function remains locked. To make changes, unlock by pressing  until  appears.

4.3 Start screen

The start screen appears for a few seconds when the unit is switched on.

4.4 Preset therapy mode

The preset therapy mode appears after the start screen; the therapy is stopped.

4.5 Changing settings

The screen of the **Suprasorb® CNP P3** is a touchscreen. The following changes can be made here by simply touching the symbols in the two toolbars (see red marking in image of screens at end of Section 4.5):

Start of therapy


The therapy mode must only be set by qualified healthcare professionals!


General settings

Locking function


4.6 Setting the therapy mode

Selecting  takes you to the following screen:

You can cancel the action by pressing  or confirm your selection by holding  down for 3 seconds to proceed to the therapy settings menu.

Choose one of the two therapy modes by selecting the corresponding view – Continuous or Intermittent – using  . The values are set in intervals of 5 for the pressure and 30 seconds for the time.


“Continuous” or “Intermittent”

In “Continuous” mode you set a constant negative pressure for the therapy.

In “Intermittent” mode, you select a negative pressure which alternates between two values. In this case, you need to specify a maximum and minimum negative pressure and set the duration of the respective setting in minutes in the time bar (under the pressure setting bar respectively).


Set the desired strength of the negative pressure in the line with the pressure symbol . In addition, set the time for which the selected pressure should be retained in the line with the clock symbol  using the   keys.

4.7 Setting the sensitivity

 <p>WARNING</p>	<p>In case of low/minimal sensitivity</p> <p>The sensitivity setting is important for ensuring that the message system is tailored to the patient’s personal wound situation. The sensitivity of the Suprasorb® CNP P3 can be set to “High”, “Medium”, “Low” or “Minimal”.</p> <ul style="list-style-type: none">• However, the sensitivity must only be set by a qualified and authorised specialist.• If the sensitivity is lowered considerably, the risk of a blockage or leak in the wound and tube system going undetected and no alarm being triggered cannot be excluded.• As such, if the sensitivity is lowered considerably, particularly careful monitoring of the therapy must be ensured.
 <p>CAUTION Concerns qualified healthcare professionals</p>	<p>Setting the therapy parameters</p> <ul style="list-style-type: none">• Before using the Suprasorb® CNP P3 on a new patient, first check that the therapy parameters set are suitable for the patient.

Not only the pressure and time but also the sensitivity need to be set. You can choose between “High”, “Medium”, “Low” and “Minimal” for this parameter on the “Therapy” screen.

The sensitivity determines how sensitively the therapy unit reacts to pressure fluctuations and how fast a message requiring attention appears on the screen. In the “High” setting, this message appears particularly quickly; in the case of “Minimal”, the delay is considerably longer.

4.8 Confirming the therapy setting

Once the setting reflects the desired therapy, confirm with .

You are automatically returned to the therapy screen, where you can now see the selected therapy settings. If you do not wish to make any changes to the therapy settings, you can return with .


In this example, the “Continuous” therapy mode is selected, the negative pressure set to -80 mmHg and the sensitivity set to “Medium”.

4.9 Starting the therapy and seal level

 <p>WARNING</p>	<p>Regular patient monitoring</p> <ul style="list-style-type: none"> Careful and regular patient monitoring is absolutely essential, particularly in patients with heavily exuding wounds, wounds with viscous or heavily sedimenting exudate and at risk of severe bleeding (see precautionary measures in Section 1.4.3). In this situation, special attention should be given to exudate accumulated in the wound dressing and the correct positioning of the wound dressing. The monitoring interval must be determined based on the wound situation, the patient's general condition and the sensitivity settings.
 <p>WARNING</p>	<p>If bleeding occurs</p> <p>The Suprasorb® CNP P3 was designed to generate negative pressure in the wound drainage system and to suction exudate.</p> <ul style="list-style-type: none"> The Suprasorb® CNP P3 was not designed to detect heavy bleeding or pooling of blood in the wound. For this reason, you should not rely on the Suprasorb® CNP P3's alarm system alone. The therapy must be regularly monitored and the exudate regularly checked for red discoloration so as to identify any risk of severe deterioration in the patient's condition. Should you observe bright red blood in the tubes or exudate pouch, the therapy must be discontinued immediately. Leave the wound dressing in place at first, stem the bleeding and seek medical assistance.
 <p>WARNING</p>	<p>Interruption of therapy</p> <ul style="list-style-type: none"> Should it prove necessary to interrupt the therapy and switch the Suprasorb® CNP P3 to Pause mode, this interruption should not be for longer than two hours at a time within a period of 24 hours. The therapy must be continued after max. two hours. If the therapy is interrupted for longer than two hours, the dressing must be changed before it can be restarted.

Press  now to start the therapy.

The screen goes green and the current pressure is shown. The seal level can be seen on the right at the bottom of the screen. It serves as a guide for how tightly the dressing is sealed.

Development of seal level:

Start phase:

Negative pressure builds up. The LED bar at the top of the device is not yet activated.

Therapy phase:

Following the start phase, one of the following screens shows that the dressing is sealed and the therapy functioning. The LED bar is green.

PLEASE NOTE

It is not possible to change the therapy settings whilst the therapy is running. You need to interrupt the therapy with  to change the settings.

Unlocking the screen

Please note that the screen lock may activate during operation. If the screen goes dark, you need only touch it to make it visible again. If you want to make changes, press  until  appears to unlock the screen.


4.10 Adjusting the language, volume, brightness, date and time

Press the  key to access the general settings menu. Here you can set the language, volume of the error messages, brightness of the screen and the date and time.

PLEASE NOTE

The “Therapy” section is exclusively reserved for qualified healthcare professionals!


Pressing the respective symbol takes you to the required settings options. If you want to save a setting, confirm with . If you want to go back without saving, press .


Pressing the “Brightness” key displays a bar showing the current brightness. This can be adjusted using the keys.

Language

To select the required language, please press “Language”. Now select the language of your choice and confirm your selection with .

Date and time

Press the “Date/time” key to access the corresponding menu. The keys can be used to set the value which is shown darker respectively.

Volume


To adjust the volume, select “Volume” and choose the volume you require with . Adjust the volume to suit the ambient conditions and ensure you hear the acoustic alarms.

<p>PLEASE NOTE</p>	<p>When the Suprasorb® CNP P3 is switched off, all the settings made after the first start are retained. The next time the device is switched on, the last mode used is displayed.</p>
---------------------------	---

Brightness

<p>PLEASE NOTE</p>	<p>Legibility of the screen</p> <ul style="list-style-type: none"> Unwanted light falling on the screen can affect its visibility. Try improving the legibility of the screen by changing the “Brightness” setting.
---------------------------	---

4.11 Viewing the chronology of the changed settings

Press the key to access the general settings menu. Use the arrow key to scroll down until “Log” appears. Confirm with and select “View”. A chronological view of the settings applied on the device is now displayed.

5 Operation by patient

PLEASE NOTE

Avoiding errors when using the device in patients' homes

- Store the **Suprasorb® CNP P3** incl. accessories in the Suprasorb® CNP Case.
- Keep children and pets away from the **Suprasorb® CNP P3** to ensure that the device is not knocked over or may fall down.
- Before connecting the power supply, ensure that the voltage of the device is the same as that of the domestic mains power supply.
- Do not use the **Suprasorb® CNP P3** in rooms with high humidity, bathrooms or shower rooms. Avoid wetting the power supply, plug and screen panel. Never immerse the **Suprasorb® CNP P3** in water or any other liquids (even when switched off).
- **Never shower or bathe with the Suprasorb® CNP P3 connected!**

5.1 Switching the device on and off

To switch the device on, press  . This key can be found below the screen on the **Suprasorb® CNP P3** (see image). To switch the device off, hold the same key down until the screen goes black.

Power supply

On/Off

Signal off

The LED of the on/off key is green when the device is switched on.

The control LED for the power supply is green when the device is connected to the mains power supply.

5.2 Screen saver and screen lock


If the screen goes dark during therapy, it can be made visible again by simply touching it. However, the editing function remains locked. To make changes, unlock by pressing  until  appears.

5.3 Start screen

The start screen appears for a few seconds when the unit is switched on.

5.4 Preset therapy mode

The screen appears with the preset parameters for your therapy. Here is an example. If the screen is orange, the therapy is paused.

5.5 Starting the therapy and switching to Pause mode

WARNING

If bleeding occurs

The **Suprasorb® CNP P3** was designed to generate negative pressure in the wound drainage system and to suction exudate.


- The **Suprasorb® CNP P3** was not designed to detect heavy bleeding or pooling of blood in the wound.
- For this reason, you should not rely on the **Suprasorb® CNP P3**'s alarm system alone. The therapy must be regularly monitored and the exudate regularly checked for red discoloration so as to identify any risk of severe deterioration in the patient's condition.
- Should you observe bright red blood in the tubes or exudate pouch, the therapy must be discontinued immediately. Leave the wound dressing in place at first, stem the bleeding and seek medical assistance.

WARNING

Interruption of therapy

- Should it prove necessary to interrupt the therapy and switch the **Suprasorb® CNP P3** to Pause mode, this interruption should not be **for longer than two hours at a time** within a period of 24 hours. The therapy must be continued after max. two hours. If the therapy is interrupted for longer than two hours, the dressing must be changed before it can be restarted.

Press the  key to start the therapy. If the screen is green, the therapy has started.

To switch the therapy to Pause mode, press . If the therapy runs for an extended period of time, the screen is

locked. To unlock it, press  until  appears. You can then switch the therapy to Pause mode using .

5.6 Adjusting the language, volume, brightness, date and time


CAUTION

Setting the therapy parameters

The therapy settings must only be changed by qualified healthcare professionals and when prescribed by the doctor.


Press the  key (see bottom line of display area) to access the general settings menu. Pressing the respective symbol takes you to the required settings options. If you want to save a setting, confirm with .

Here you can set the language, volume of the error messages, brightness of the screen and the date and time. The therapy settings must only be changed by qualified healthcare professionals!

If you want to go back without saving, press .

Language

To select the required language, please press "Language". Now select the language of your choice and confirm your selection with .

Volume

To adjust the volume, select “Volume” and choose the volume you require with **-** **+**. Adjust the volume to suit the ambient conditions and ensure you hear the acoustic alarms.


Brightness

 <p>PLEASE NOTE</p>	<p>Legibility of the screen</p> <ul style="list-style-type: none">Unwanted light falling on the screen can affect its visibility. Try improving the legibility of the screen by changing the “Brightness” setting.
--	---

Pressing the “Brightness” key displays a bar showing the current brightness. This can be adjusted using the **-** **+** keys.

Date and time

Press the “Date/time” key to access the corresponding menu.

Press once to select the day/month/year or time and then use the **-** **+** keys to adjust the respective value.

5.7 If the doctor has granted permission to shower

 <p>WARNING</p>	<p>Interruption of therapy</p> <ul style="list-style-type: none">Should it prove necessary to interrupt the therapy and switch the Suprasorb® CNP P3 to Pause mode, this interruption should not be for longer than two hours at a time within a period of 24 hours. The therapy must be continued after max. two hours. If the therapy is interrupted for longer than two hours, the dressing must be changed before it can be restarted.
---	---

If the doctor has granted you permission to shower, proceed as follows:

1. Close the tube clamps on the exudate pouch tube and drainage pad tube.
2. Switch the therapy to Pause mode (see Section 5.5).
3. Disconnect the two tubes and leave the device next to your bed whilst you shower.
4. After showering, connect the two tubes again (see Section 3.3) and restart the therapy (see Section 5.5).

6 Handling error messages

6.1 Function test

Proceed as described in Section 3.2.

6.2 Device does not start

Error	Possible cause	Remedy
Device does not start	<ul style="list-style-type: none">• Battery is empty• Tube clamp is closed• Overflow prevention system is blocked (exudate pouch is full)• Internal filter blocked• Device is still in <i>Settings</i> mode	<ul style="list-style-type: none">• Connect power supply• Check tubes are connected correctly• Replace exudate pouch• Please contact the Service department• Confirm selection (see Section 3.2) and start the device

PLEASE NOTE

Contact L&R or the service partner authorised by L&R if an error occurs that you cannot resolve with the measures described here.

You can reach the L&R CNP service hotline 24 hours a day in the following countries:

Germany: +49 (0)800 6647435
Austria: +43 (0)800 808106
Switzerland: +41 (0)848 343434

6.3 Error messages on the screen and solutions

When there is an error message, the LED bar goes orange, the screen background goes red and an acoustic alarm sounds. The volume can be changed in the settings (see Section 4 or 5 respectively). In addition, the alarm tone can be muted for two minutes using the  key on the housing.

PLEASE NOTE

Please proceed as follows:

1. When an error message appears, mute the alarm tone
2. Resolve the problem
3. Then return with ; the therapy continues automatically

Error message	Status	Possible cause	Remedy
	<p>A signal sounds. Therapy is interrupted.</p>	<p>Container is not connected. Connections open.</p>	<p>The container and exudate pouch may not be connected to the device properly. The cause may also be a collapsed exudate pouch. Or a connection between the tubes may be open. Please check the points above and re-establish all connections if necessary.</p>
	<p>A signal sounds. Therapy continues to run in the background.</p>	<p>Wound dressing or connections not tightly sealed.</p>	<p>The tubes might not be connected properly or the dressing itself is not sealed tight. Please check both possibilities, reconnect the tubes using the single-piece connector or improve the seal on the dressing.</p>
	<p>A signal sounds. Therapy continues to run in the background.</p>	<p>Exudate flow interrupted, tube blocked, internal filter blocked.</p>	<p>A tube clamp may have accidentally been closed or there may be a blockage in the tubes. Please open all the tube clamps or try to remove the blockage in the tubes. If necessary, change the full exudate pouch. However, if the internal charcoal/bacterial filter is blocked, the device needs to be reprocessed properly by L&R or an authorised L&R service partner!</p>
	<p>A signal sounds. Therapy continues to run in the background.</p>		<p>If the message “Exudate pouch is full” appears, it is high time to insert a new exudate pouch. Replace the pouch in accordance with the instructions for use (see Section 3.4).</p>

Error message

Status

Possible cause

Remedy

 <p>Battery low</p> <p>Please connect the device to the power supply.</p>	<p>A signal sounds. Therapy continues to run in the background.</p>	<p>Low battery level.</p>	<p>If the battery level is low, please connect the device to the mains power supply using the supplied power supply.</p>
 <p>Battery empty</p> <p>Please connect the device to the power supply immediately.</p>	<p>A signal sounds. Therapy is interrupted.</p>		<p>If the battery level is empty, please connect the device to the mains power supply immediately using the supplied power supply.</p>
 <p>No therapy</p> <p>Please start the therapy or turn the device off.</p>	<p>A signal sounds.</p>	<p>Therapy has not been started.</p>	<p>You now have the option of starting the therapy or switching the device off.</p>
 <p>Internal error</p> <p>Please turn the device OFF and ON again. If the error is still displayed, please contact your CNP specialist.</p>	<p>A signal sounds. Therapy is interrupted.</p>		<p>A simple trick often helps in such cases: Switch the device off and on again, then check whether it is functioning normally. If this does not help, please contact your CNP specialist or the L&R hotline (see Section 11).</p>

PLEASE NOTE

If you cannot resolve an error message successfully with the measures outlined above, contact the L&R hotline (see Section 11) or your CNP specialist.

7 Cleaning of Suprasorb® CNP P3 therapy unit

CAUTION

Avoiding infections

Infectious and pathogenic microbes in the exudate can pose health risks.

- Wear appropriate disposable gloves when **replacing the container**.
- Replace the **exudate pouch** including the exudate pouch tube in accordance with the respectively applicable hygiene regulations, at the latest every week, before a change of patient or when full.
- **Reprocessing is required whenever there is a change of patient.** Reprocessing means that all the single-use products are disposed of and the **Suprasorb® CNP P3** therapy unit wipe-down disinfected. The **Suprasorb® CNP P3** therapy unit must only be opened by L&R or an authorised L&R service partner.
- The exudate and all components contaminated with it must be disposed of **in accordance with the facility's policies**.

CAUTION

Cleaning the device

- Switch off the device prior to cleaning/disinfection by disconnecting the power supply from the mains power and then from the power socket on the **Suprasorb® CNP P3**.
- Disconnection from the power supply is only guaranteed when the charging cable is removed from the charging port.

PLEASE NOTE

Use of disinfectants

- We recommend wearing appropriate disposable gloves when performing disinfection.
- Please observe the disinfectant manufacturer's specifications.

Disinfection during the therapy is not expressly required when using the device on a single patient at home.

Disinfection is generally required when the device is used in an inpatient setting.

7.1 Cleaning and disinfection of the device surfaces and the container

- Clean the device regularly and disinfect the surfaces at least once a week.
- All devices can be wiped with a damp, lint-free cloth.

When the device has been in use for an extended period, the plastic components of the housing may exhibit slight discolouration. This does not have any effect on the functioning of the device.

If liquids come into contact with the inside of the device, the unit must be checked by L&R or an L&R authorised service partner.

7.2 Maintenance and service

 <p>CAUTION</p>	<p>Use of the carry bag</p> <ul style="list-style-type: none">• The Suprasorb® CNP Carry Bag for mobile use supplied with the Suprasorb® CNP P3 is intended for single use per patient. It should not be used for other patients but be disposed of once used.
 <p>PLEASE NOTE</p>	<p>Cleaning of the therapy unit</p> <ul style="list-style-type: none">• Aggressive cleaning agents may damage the Suprasorb® CNP P3 and its accessories. Do not use any disinfectants containing acetone.• Follow the manufacturer's instructions for use for the disinfectants used, especially with regard to the compatibility with the materials and surfaces as well as the concentration specifications.• Follow the cleaning and care instructions (Section 7). L&R recommends Incidin LIQUID for wipe disinfection of the device.
 <p>PLEASE NOTE</p>	<p>The Suprasorb® CNP P3 therapy unit must undergo a check in accordance with DIN EN 62353 performed by L&R or an L&R authorised service partner every 18 months at the latest. L&R offers customers and partners rapid and professional reprocessing as well as the performance of necessary tests and checks.</p>

The **Suprasorb® CNP P3** therapy unit requires no maintenance when used in accordance with the instructions for use.

Always perform a visual inspection and function test prior to use. Do the same with the device's accessories.

The device may only be opened and repaired by L&R or an authorised specialist in accordance with the service documentation provided by the manufacturer and taking appropriate technical and hygienic safety precautions.

The device may be returned to L&R directly or via your CNP specialist.

All accessories must be cleaned and disinfected prior to return. The device itself must be treated with a surface disinfectant.

L&R does not guarantee the error-free functioning of the **Suprasorb® CNP P3** therapy unit nor does it accept liability for any material damage or personal injuries sustained if:

- accessories and spare parts are used which are not original;
- the information in these instructions for use is not heeded;
- assembly, setup, modifications, updates or repairs are not performed by L&R or an authorised specialist;
- the security seal is removed or broken.

8 Transport, storage and disposal

Please read “Important information on your safety”, Section 1.5.

8.1 Decontamination before shipping

 <p>CAUTION</p>	<p>Avoiding infections</p> <p>Infectious and pathogenic microbes in the exudate can pose health risks.</p> <ul style="list-style-type: none">• Wear appropriate disposable gloves when replacing the container.• Replace the exudate pouch including the exudate pouch tube in accordance with the respectively applicable hygiene regulations, at the latest every week, before a change of patient or when full.• Reprocessing is required whenever there is a change of patient. Reprocessing means that all the single-use products are disposed of and the Suprasorb® CNP P3 therapy unit wipe-down disinfected. The Suprasorb® CNP P3 therapy unit must only be opened by L&R or an authorised L&R service partner.• The exudate and all components contaminated with it must be disposed of in accordance with the facility's policies.
 <p>CAUTION</p>	<p>Use of the carry bag</p> <ul style="list-style-type: none">• The Suprasorb® CNP Carry Bag for mobile use supplied with the Suprasorb® CNP P3 is intended for single use per patient. It should not be used for other patients but be disposed of once used.
 <p>CAUTION</p>	<p>Shipping instructions</p> <ul style="list-style-type: none">• The therapy unit must be sent in the supplied case, which must also be protected by an additional box.• Used exudate pouches must not be sent with the device but should be disposed of properly.

Before the **Suprasorb® CNP P3** therapy unit is passed on to the next user, it must be professionally reprocessed by L&R or an authorised specialist.

The reprocessing must be performed in accordance with the Medical Devices Operator Ordinance (MPBetreibV), the Medical Devices Act (MPG) and the manufacturer's specifications. L&R offers its customers and partners rapid and professional reprocessing as well as the performance of necessary tests and checks.

The **Suprasorb® CNP P3** therapy unit must be cleaned and disinfected before dispatch to the L&R Service Centre. Please refer to Section 7.

8.2 Storage

- Store the **Suprasorb® CNP P3** therapy unit in accordance with the specifications in Technical Data (Section 9).
- Charge the **Suprasorb® CNP P3** therapy unit's battery before storing the device. This ensures that it can be used at any time.
- If the **Suprasorb® CNP P3** therapy unit is not used for an extended period of time (several months), the battery should be fully recharged!

8.3 Disposal

- At the end of the product's service life, the components of the device must be disposed of in accordance with the applicable regulations.
- Ensure that all parts are clean and separate materials accordingly for recycling.
- The housing material is labelled with a material symbol and can be completely recycled.
- Decontaminate the device and container prior to disposal. For information on disposal of the accessories (exudate pouch, drainage pad, deep drain, Y connector), please refer to the instructions for use supplied with the respective products.
- In accordance with EU Directive 2002/96/EC on Waste Electrical and Electronic Equipment Directive (WEEE) and 2002/95/EC on the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment (RoHS II), the device must not be disposed of as household waste.
- The device and accessories can be disposed of via L&R or an authorised service partner.
- Please observe any applicable national disposal regulations.

9 Technical data of the Suprasorb® CNP P3 therapy unit and EMC information

Aggregate suction power	11 l/min
Negative pressure	-20 to -200 mmHg*; conversion factor: 1 kPa≈7.5 mmHg
Container	Suprasorb® CNP Container 250 ml Suprasorb® CNP Container 1,000 ml
Power supply rated voltage	100-240 V AC primary / 15 V DC secondary
Maximum load current	1 A
Power supply line frequency	50/60 Hz
Power consumption	15 W (charging and operation)
Battery, rechargeable	10.8 V, 2.7 Ah – lithium ion
Charging time when battery is empty	6-7 hours
Charging time with approx. 50% residual charge	3-3.5 hours
Dimensions (H x W x D)	175x155x70 mm with Suprasorb CNP Container 250 ml 175x260x70 mm with Suprasorb CNP Container 1,000 ml
Weight (therapy unit)	1.15 kg with Suprasorb CNP Container 250 ml 1.38 kg with Suprasorb CNP Container 1,000 ml
Operating time	Mains: Continuous operation Battery: approx. 10-20 hours
Degree of protection in accordance with EN 60601-1	Type BF 
Risk class as per 93/42/EEC, Annex IX, Rule 11	Ila
Protection class in accordance with EN 60601-1	II
CE marking	CE0123
Noise emissions	40 dB (A)
Ambient conditions	Transport/storage: -10°C to +60°C Operation: +5°C to +40°C Recommended charging temperature: +15 C to +30 C 5 to 80% relative humidity, non-condensing Atmospheric pressure: 860 hPa to 1,060 hPa
REF	31780 34780

* Maximum negative pressure fluctuations:

The maximum value refers to sea level.

It may vary depending on altitude and atmospheric pressure. The multiplication factors listed in the table are employed to set the maximum negative pressure based on the altitude and atmospheric pressure.

Location	0 metres above sea level	500 metres above sea level	1,000 metres above sea level	2,000 metres above sea level
Multiplication factor	1.0	0.94	0.88	0.78

Table 1

Guidelines and manufacturer's declaration – Electromagnetic emissions

The **Suprasorb® CNP P3** therapy unit is intended for operation in an ELECTROMAGNETIC ENVIRONMENT as described below.

The customer or user of the **Suprasorb® CNP P3** therapy unit must ensure that it is operated in such an environment.

Measurement of emitted interference

Conformity

Electromagnetic environment – guidance

HF emissions in accordance with CISPR 11	Group 1	<p>The Suprasorb® CNP P3 therapy unit only uses HF energy for its internal workings. As such, its HF emissions are very low and it is unlikely that adjacent electronic devices will be affected.</p> <p>The Suprasorb® CNP P3 therapy unit is intended for use in all facilities including residential buildings and those directly connected to the mains power grid which also supplies residential buildings.</p>
HF emissions in accordance with CISPR 11	Class B	
Harmonic current emissions in accordance with IEC 61000-3-2	Class A	
Emission of voltage fluctuations/flickers in accordance with IEC 61000-3-3	Conforms	

Table 2

Guidelines and manufacturer's declaration – Electromagnetic immunity

The **Suprasorb® CNP P3** therapy unit is intended for operation in an electromagnetic environment as described below.

The customer or user of the **Suprasorb® CNP P3** therapy unit must ensure that it is operated in such an environment.

Immunity tests	IEC 60601 test level	Conformity level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) in accordance with IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be made of wood or concrete or covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be min. 30%.
Electrical fast transient/bursts immunity in accordance with IEC 61000-4-4	± 2 kV for power cables ± 1kV for input and output cables	± 2 kV for power cables ± 1kV for input and output cables	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Surges in accordance with IEC 61000-4-5	± 1 kV voltage phase-phase ± 2 kV phase-earth	± 1 kV voltage phase-phase ± 2 kV phase-earth	The quality of the supply voltage should correspond to that of a typical business or hospital environment.
Voltage dips, short interruptions and voltage variations in accordance with IEC 61000-4-11	<5% UT (>95% dip of the UT) for 0.5 cycle 40% UT (60% dip of the UT) for 5 cycles 70% UT (30% dip of the UT) for 25 cycles <5% UT (>95% dip of the UT) for 5 seconds	<5% UT (>95% dip of the UT) for 0.5 cycle 40% UT (60% dip of the UT) for 5 cycles 70% UT (30% dip of the UT) for 25 cycles <5% UT (>95% dip of the UT) for 5 seconds	The quality of the supply voltage should correspond to that of a typical business or hospital environment. If the user of the Suprasorb® CNP P3 therapy unit requires it to continue functioning even if the power supply is interrupted, we recommend connecting the device to a power supply without interruptions or running it off the battery.
Magnetic fields at the supply frequency (50/60 Hz) in accordance with IEC 61000-4-8	3 A/m	3 A/m	Magnetic fields at the supply frequency should be equivalent to the typical values encountered in business and hospital environments.
NOTE: UT is the mains alternating voltage prior to application of the test level.			

Table 3

Guidelines and manufacturer’s declaration – Electromagnetic immunity

The **Suprasorb® CNP P3** therapy unit is intended for operation in an electromagnetic environment as described below.

The customer or user of the **Suprasorb® CNP P3** therapy unit must ensure that it is operated in such an environment.

Immunity tests	IEC 60601 test level	Conformity level	Electromagnetic environment – guidance
Conducted HF disturbances in accordance with IEC 61000-4-6	3 V	3 V	<p>Portable and mobile wireless devices should not be kept closer to the Suprasorb® CNP P3 therapy units, including all cables, than the recommended protective distance calculated using the applicable equation for the transmission frequency.</p> <p>Recommended protective distance:</p> $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where P is the maximum nominal capacity of the transmitter in watts (W) in accordance with the transmitter manufacturer’s specifications and d is the recommended protective distance in metres (m). The field strength of stationary transmitters should be less than the conformity level at all frequencies according to a non-site^a study.^b Interference is possible in the vicinity of devices bearing the following symbol.</p> 
Radiated HF disturbances in accordance with IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1: The higher frequency range applies at 80 MHz and 800 MHz.

NOTE 2: These guidelines may not be applicable in all cases. The transmission of electromagnetic parameters is influenced by the extent to which they are absorbed and reflected by buildings, objects and people.

^a In theory, it is not possible to determine the field strength of stationary transmitters, such as base stations of radio phones and land mobile radio services, CB radio stations, AM and FM radio and TV stations in advance with accuracy. A site study of the prevalent electromagnetic phenomena should be taken into consideration for the determination of the electromagnetic environment with regard to stationary transmitters. If the field strength measured at the location at which the **Suprasorb® CNP P3** therapy unit is used exceeds the above conformity levels, the device should be monitored to verify that it functions as intended. If unusual performance characteristics are observed, additional measures such as a modified alignment or a different location for the device may prove necessary.

^b The field strength should be less than 3 V/m over the frequency range from 150 kHz to 80 MHz.

Table 4

Recommended safe distances between portable or mobile HF telecommunication devices and the **Suprasorb® CNP P3** therapy unit

The **Suprasorb® CNP P3** therapy unit is intended for operation in an electromagnetic environment as described below in which HF interference is controlled.

The customer or user of the **Suprasorb® CNP P3** therapy unit can help avoid electromagnetic interference by ensuring maintenance of the minimum distance between portable and mobile HF telecommunication devices (transmitters) and the **Suprasorb® CNP P3** – depending on the output power of the communication device, as specified below.

Nominal capacity of transmitter	Protective distance depending on transmitter frequency in m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1}\right]\sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33

For transmitters with maximum nominal power not given in the table above, the recommended protective distance (d) in metres (m) can be calculated using the formula in the respective column, where P is the transmitter's maximum nominal power in watts (W) in accordance with the transmitter manufacturer's specifications.

NOTE 1: The higher frequency range applies at 80 MHz and 800 MHz.

NOTE 2: These guidelines may not be applicable in all cases. The transmission of electromagnetic parameters is influenced by the extent to which they are absorbed and reflected by buildings, objects and people.

10 Ordering options

Ensure that you always have sufficient supplies of accessories and spare parts available.

L&R only guarantees the operational security of the device if the specified original accessories and spare parts are used.

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www.lohmann-rauscher.com

Hotlines:

The following countries operate an L&R CNP service hotline 24/7.

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