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Instructions for Use

Product description

Laparotomy Swabs RK, non-sterile, consist of 17-thread count cotton gauze in accordance with EN 14079. The included X-ray detectable strip is made out of polyvinyl chloride (PVC) and contains barium sulphate. The laparotomy swabs have a blue loop, which is also made out of cotton.

Product composition

100% cotton, X-ray detectable strip (blue, polyvinyl chloride, barium sulphate)

Intended purpose

Laparotomy Swabs made of cotton gauze are used in sterile condition for grasping and retaining organs and body tissues as well as for haemostasis and absorbing fluids during surgical interventions.

Usage information/Precautions

This product is non-sterile. Laparotomy Swabs must be sterilised prior to surgically invasive use; resterilisation is not permitted. The product is intended for single use only.

Sterilisation instructions

1. Device requirements

Only sterilisers, which comply to the requirements of the norms EN 285 or EN 13060 and which use a pre-vacuum system for the sterilisation are allowed to be used.

The use of the gravitational method is not allowed in the EU.

2. Sterilisation preparation

Before the product is sterilised, a suitable sterile barrier system must be selected, validated and used by the person responsible for the sterilisation (sterilising technician) in accordance with ISO 11607. The sterile barrier must be big enough for the product so that the seal is not under tension. The retail packaging is unsuitable for sterilisation.

3. Sterilisation

The sterilisation procedure for the product must be validated by means of moist heat at 134°C for 7 minutes in accordance with ISO 17665. Do not use after sterilisation if package is damaged or accidentally opened.

4. Storage

After removing the required quantity of Laparotomy Swabs, non-sterile, seal the packaging and store dry and protected from dust. Sterilised products must be stored dry and protected from dust in microbe-resistant sterile packaging.

Additional information

The above instructions were assessed by Lohmann & Rauscher as suitable for preparation for sterilisation. The sterilising technician is responsible for ensuring that the sterilisation carried out with the employed equipment, materials and staff achieves the required results. To this end, validation of the sterilisation process and routine monitoring of the procedure is normally required. Deviations from the procedure are the responsibility of the sterilising technician. National requirements for steam sterilisation must be considered by the user.

Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

General instructions

In case of a serious incident, contact the manufacturer and the responsible health authority.

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