Suprasorb® CNP

Y-Konnektor
Connecteur en Y
Y-Connector
Conector en Y
Conetor em Y
Connettore a Y
Y-connector
Y-konnektor

Y-koppling
Y-konektor
Y-konnektor
Y-οβραζεν konektor
Triguba jungtis
Y konnektörü
Yłacze Y
Złącze Y
Y-konektor
Y konektor

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The instructions for use for the Suprasorb® CNP P3 system must also be complied with, in particular with regard to indications, contraindications, side effects, warnings and precautions.

Product description and features
The Suprasorb® CNP Y-Connector consists of the connector body, three sliding sleeves, the sealing sleeves and two expanding seals.
The Suprasorb® CNP Y-Connector is connected to the Suprasorb® CNP Exudate Pouch by the main line. The two remaining lines can be connected to two Suprasorb® CNP Deep Drains, two Suprasorb® CNP Drainage Pads or a Suprasorb® CNP Deep Drain and a Suprasorb® CNP Drainage Pad.

Product composition
Connector body: Polypropylene
Sliding sleeve: Acrylonitrile, butadiene, styrene
Sealing sleeve: Polypropylene, thermoplastic elastomers
Expanding seal: Polypropylene, cellulose, superabsorbent polymer, Hot Melt

Intended purpose
The Suprasorb® CNP Y-Connector is used in negative pressure therapy for the simultaneous connection of two lines from the wound to the Suprasorb® CNP P3 Therapy Unit.

Indications
Not applicable.

Contraindications
Known allergy and / or hypersensitivity to any of the product components.

Side effects
In very rare cases skin irritations and / or allergies may occur.

Warnings and Precautions
• Monitor the therapy closely when employing the Suprasorb® CNP Y-Connector. • The signals emitted by the connected Suprasorb® CNP P3 Therapy Unit are not always clear when the Suprasorb® CNP Y-Connector is employed.
• Do not use the Suprasorb® CNP Y-Connector on wounds of differing aetiology, as such actions may pose the risk of cross-contamination.
• The Suprasorb® CNP Y-Connector must only be used in combination with the Suprasorb® CNP P3 Therapy Unit.
• This product is intended for use exclusively by healthcare professionals.
• El conector en Y Suprasorb® CNP solo debe utilizarse en combinación con la unidad de terapia Suprasorb® CNP P3.
• El producto debe ser utilizado exclusivamente por personal médico especializado.

Application instructions

1. Application
a. Prior to use, visually inspect the sterile barrier system for breaches to the packaging integrity.
b. Open the peel pouch at the place provided for this purpose and remove the product, keeping it sterile.
c. Open all three ends of the Suprasorb® CNP Y-Connector by pushing the sliding sleeves on the Suprasorb® CNP Y-Connector TOWARDS EACH OTHER in the direction indicated by the arrows. (See Figure 1)
d. Ensure that any existing connectors are removed from the ends of the tubes that you intend to connect to the Suprasorb® CNP Y-Connector. (See Figure 2)
e. Insert the exudate pouch tube into the longer branch of the Suprasorb® CNP Y-Connector and the tubes from the wound into the two shorter branches.
f. Insert the tubes TWO CENTIMETRES into the connector. Push the connector in Y Suprasorb® CNP (consulte la figure 2).
g. Close all three ends of the Suprasorb® CNP Y-Connector by pushing the sliding sleeves on the Suprasorb® CNP Y-Connector APART in the direction indicated by the arrows. (See Figure 4)

2. Change
The Suprasorb® CNP Y-Connector is changed at each dressing change.

3. Disposal
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems. Additionally, please observe any applicable national disposal regulations.

General instructions
In case of a serious incident, contact the manufacturer and the responsible health authorities.

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