Idealbinde

Idealbinde steril

Kurzzügige Idealbinde nach DIN 61 632
Bande à allongement courte Idealbinde selon la norme DIN 61 632
Short-stretch Idealbinde in accordance with DIN 61 632
Venda de tracção corta Idealbinde según DIN 61 632
Idealbinde de baixa elasticidade conforme a norma DIN 61 632
Idealbinde ad allungamento limitato conforme alla norma DIN 61 632
Idealbinde met korte rek volgens DIN 61 632
Kortstræks-Idealbinde i henhold til DIN 61 632
Lågelastisk Idealbinde enligt DIN 61 632
Idealbinde elastické obinadlo s krátkým tahem dle DIN 61 632
Ovinadlo Idealbinde s krátkým ťahom podľa normy DIN 61 632
Idealbinde rövid megnyúlású pólya a DIN 61 632 szerint
Bandaž o krótkim naciagu Idealbinde zgodny z DIN 61 632
DIN 61 632'ye uygun kısa çekişli Idealbinde

Lenkideal®

Kurzzugbinde
Bande à allongement court
Short Stretch Bandage
Venda de tracção corta
Ligadura de baixa elasticidade
Benda ad allungamento limitato
Korte-rekzwachtel
Kortstræksbind
Lågelastisk binda
Obinadlo s krátkým tahem
Ovinadlo s krátkým ťahom
Rövid megnyúlású pólya
Bandaž o krótkim naciagu
Kısa çekişli bandaj

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en Instructions for Use

Product description and performance characteristics:
• Idealbinde: Stretch approx. 90%
• Lenkideal®: Stretch approx. 85%

Product composition:
• Idealbinde: 100% cotton
• Lenkideal®: 64% cotton, 34% polyamide, 2% elastane

Intended purpose:
For fixation, support and relief, for immobilisation of limbs and for compression.

Indications:
phlebological and lymphological indications: varicosis, deep venous insufficiency, thrombophlebitis, phlebothrombosis, post-thrombotic syndrome, Widmer's stage I to III b chronic venous insufficiency, venous and mixed ulcers, oedema (phleboedema, lymphoedema, lipoedema, post-traumatic, post-operative, during pregnancy), after sclerotherapy and venous surgery.
contusions, sprains, tendon injuries and joint dislocations

Contraindications:
Advanced peripheral arterial occlusive diseases, congestive heart failure, septic phlebitis, phlegmasia cerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to any of the product components.

Side effects:
In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions:
• Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
• Slight bluish tinge to the toes is normal, this should disappear as soon as the patient starts to walk.
• If the symptoms do not disappear (pain and bluish tinge), the bandage must be removed and applied again with less compression.
• The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage.
• For use with ABPIs over 0.8. When ABPI is between 0.5 and 0.8, apply under medical supervision only.
• Use adhesive strips to fix the bandage in place. Do not fix in place with bandage clips because these can cause injury.
• If appropriate, an underlying bandage is advisable – depending on the condition of the skin.
• Known allergy and / or hypersensitivity to any of the product components.

General instructions:
• Before applying the bandage, clean/debride any wound that is present and cover with a wound dressing appropriate for the wound status.
• It is recommended to apply the bandage in the morning before the patient gets up or after the patient has rested the leg in a raised position for a prolonged period and to then encourage the patient to walk.
• The width of the bandage should correspond to the diameter of the limb to be bandaged.
• To achieve maximum joint movement, it is recommended to apply the bandage with the foot at an angle of 90° to the leg.
• The bandage must be applied so that the pressure from the bandage decreases from the foot to the knee and, where applicable, to the gluteal fold.
• The bandage is unrolled with the rolled-up portion on top. Keep the roll close to the leg/arm.

Usage information:
• Idealbinde and Lenkideal® can be applied using different compression techniques (illustrations A1-A11 or B1-B10).
• The bandaging technique is chosen in consultation with the treating physician.

Care recommendation:
Washing may cause the ends of the bandage to fray. It will retain its short-stretch characteristics and function. The maximum stretch of the bandage may increase with washing. Bear this in mind when applying the bandage.
Idealbinde: Can be washed up to 5 times at 95°C (standard commercially available mild detergent/delicates cycle).

Disposal:
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General instructions:
In case of a serious incident, contact the manufacturer and the responsible health authorities.

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