de **Gebruuchsanweisung**

Suprasorb® Liquacel

Hydroaktiver Wundverband / Hydroaktive Tamponade

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en Instructions for Use

Product Description

Suprasorb® Liquacel wound dressings are a soft, conformable non-woven dressings made from a blend of gel forming cellulose ethyl sulfonate (CES fibres) and strengthening cellulose fibres(s). When it comes into contact with wound exudate or blood the absorbent dressing forms a gel which creates a moist wound environment that may support autolytic debridement and the healing process. Through the gel formation, debris and any bacteria found in the wound exudate can be retained inside the fibre dressing and removed when the dressing is changed. When dry, Suprasorb® Liquacel wound dressings can easily be cut to the size of the wound. Even when the dressing is moist, as a gelled fibre, its structure remains intact. The high vertical absorption of exudate into the fibre dressing protects the wound environment and the wound edge¹, thus supporting the healing process. The product is for a single use only.

Indications for use

This product is to be used for wounds with moderate to heavy exudate. It may be used for the treatment of acute or chronic wounds;

- Leg ulcers, pressure ulcers (Stagell to IV) and diabetic ulcers
- Surgical wounds (e.g. post-operative, wounds left to heal by secondary intent and donor sites)
- Partial thickness burns
- Traumatic wounds (e.g. abrasions and lacerations)
- Exuding oncology wounds

provided they are:

 moderately to heavily exuding and

 superficial or

 deep:

The product can also be used under compression. The Suprasorb® Liquacel wound dressing has a supportive effect in that it protects the wound edge and the surrounding skin from maceration.¹

Contraindications

The product should not be used on people who are sensitive or allergic to the dressing and its components.

Application

Preparing the Wound

Clean the wound carefully. The skin surrounding the wound should be clean and dry.

Applying the Wound Dressing

- The size of the dressing should match the wound area. Keep the product away from the wound while cutting it to size.
- Place the dressing on the moist wound, overlapping the wound edges by approximate 1cm. If there is a low amount of exudate present, moisten the dressing the physiological saline solution (0.9%).
- For cavity wounds, use the rope format. Pack wounds lightly and leave approximate 2,5cm overlapping the wound edges.
- Apply a suitable moisture-retaining secondary dressing to keep the product in place.

Changing the Dressing

The dressing should be changed when medically indicated (e.g. when the dressing has reached its absorbent capacity or whenever good wound care practice dictates a change is needed). The interval between changes should be no more than 7 days. The product can be removed using e.g. sterile forceps. Should the dressing adhere to wounds with lower exudate levels, moisten it with physiological saline solution (0.9 %) before changing the dressing so that the healing process is not disturbed. Any gel residue on the wound should be removed when cleansing the wound. Wound cavities in particular should be well irrigated. The treatment should be discontinued after 30 days. Dispose in accordance with local guidance.

Caution

- If signs of infection are detected clinically the healthcare professional responsible for the treatment must decide on the next course of action.
- If signs of an allergic reaction are detected, discontinue the use of the product and seek advice from a healthcare professional on the next course of action.
- Product is sterilised by irradiation and must not be re-sterilised.
- The product remains sterile unless the package is opened or damaged.
- Suprasorb® Liquacel is not intended to be used on wounds with severe bleeding, inside internal body cavities or closed wounds.

¹ The properties mentioned above refer to the current in vitro data.

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