Suprasorb® CNP

Tiefendrainage

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Deep drain

Drenaje profundo

Drenaggio profundo

Dieptedrainage

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Dijpdrainage

Dybðdrenening

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Mély drenázs

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en Instructions for use

The instructions for use for the Suprasorb® CNP P3 system must also be complied with, in particular with regard to indications, contraindications, side effects, warnings and precautions.

Product description and performance characteristics
The Suprasorb® CNP Deep Drain is a 90-centimetre-long tube in which an X-ray contrast strip, a central exudate lumen and three ventilation lumens are integrated. Additional components of the Suprasorb® CNP Deep Drain include a tube clamp and a sleeve secured with adhesive on the end section. The supplied connector establishes the secure connection to the tubing of the Suprasorb® CNP P3 therapy system. The Suprasorb® CNP Deep Drain is intended for perforation of the wound from the wound features over a length of 30 centimetres. The remaining 60 centimetres of the silicone tube are not perforated.

Product composition
Tube: silicone
X-ray contrast strip: 60% silicone, 40% barium sulphate
Tube clamp: polypropylene
Sleeve: PVC
Connector: polypropylene, acrylonitrile butadiene styrene, thermoplastic elastomers
Adhesive: cyanoacrylate

Intended purpose
The drainage tube is placed inside a wound filler during negative pressure therapy with the Suprasorb® CNP P3 system in a wound area to transport exudate from the wound to the exudate pouch.

Indications
Not applicable.

Contraindications
Known allergy and / or hypersensitivity to any of the product components.

Side effects
In very rare cases skin irritations and / or allergies may occur.

Warnings and Precautions
• The Suprasorb® CNP Deep Drain must only be used in combination with the Suprasorb® CNP P3 therapy unit.
• This product is intended for use exclusively by healthcare professionals.

3. Dressing change:
Dressing change intervals depend on the patient and wound condition. The wound filler used and the therapy system, as well as the assessment of the treating physician.

4. Disposal:
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (AWfK) (European Waste Catalogue).

General instructions
In case of a serious incident, contact the manufacturer and the responsible health authorities.

Application instructions
When opening the sliding folding box, ensure that the Suprasorb® CNP P3 system is not damaged.

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