Curafix® i.v. Control

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en Instructions for use

Product description and features:

Curafix® i.v. Control is a retention dressing for catheters that is particularly suitable for central access. It consists of white non-woven material, adhesive material, film and cover paper. The product also includes a counterfixation element. The top part of the product has an integrated viewing window that allows monitoring of the cannula/catheter insertion site and the insertion channel. The dressing can be washed with water, soapy water and iodine- or octenidine-containing disinfectants based on alcohol or water.

Product composition:

Polyester, polyacrylate, polyurethane, glassine paper with polydimethylsiloxane coating

Intended use:

For the retention of cannulas and catheters on the skin and as a mechanical barrier to protect the puncture site.

Indications:

Not applicable.

Contraindications:

Known allergy and/or hypersensitivity to any of the product components.

Side effects:

In very rare cases skin irritations and/or allergies may occur.

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Instructions for use:

Application is performed by healthcare professionals.

1. Preparation:

After thorough cleaning and disinfection of the puncture site and the surrounding skin, the dressing is applied to the dried skin.

2. Applying the retention dressing (after puncture):

Before opening the sterile barrier, check for obvious damage.

a) Open the packaging at the point indicated and remove the product in a sterile manner (Fig. 1).

b) Remove the large cover paper, on which the counterfixation element is also attached, by pulling it downwards and initially place back in the packaging (Fig. 2).

c) Place the dressing with the V-shaped incision around the catheter or lumen exiting the skin so that the puncture site is covered by the transparent viewing window (Fig. 3).

d) Now remove the cover paper on the viewing window and mould the dressing in place (Fig. 4).

e) Then remove the large cover paper from the packaging, detach the counterfixation element and place and fix around the catheter or lumen in the opposite direction to the previously applied dressing (Fig. 5).

f) Finally, pull off the grip tab (small cover paper) and press on the counterfixation element (Fig. 6).

3. Dressing change:

The intervals at which the dressing must be changed are determined by the treating physician. These depend on the patient's individual skin condition and on the underlying medical findings.

Replace dressings which are soaked through, soiled or no longer adhere securely. Also replace the dressing if weeping pockets develop under the dressing.

4. Additional note:

To detach Curafix® i.v. Control, the use of dressing-loosening fluids (based on alcohol, saturated carbohydrates or siloxanes) is recommended.

5. Disposal:

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

6. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

7. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

8. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

9. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

10. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

11. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

12. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

13. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

14. General instructions:

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

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