Rosidal® soft

Ulcus cruris Kompressionssystem Système de compression pour ulcère de la jambe Leg ulcer compression system Sistema de compresión para úlceras de la pierna Sistema de compressão para úlcera de perna Sistema compressivo per l'ulcera vascolare Compressiesysteem voor ulcus cruris Kompressionssystem til ulcus cruris Kompressionssystem for bensår Kompresní systém pro léčbu ulcus cruris Kompresný systém na terapiu ulcus cruris Lábszárfekély kompressziós rendszer System do kompresji w przypadku owrzodzenia podudzi Kompresijski sistem za golenje razjede Kompresivni sustav za potkoljenični ulkus Bacak ülseri kompresyon sistemi

Rosidal® K

50x

Rosidal® soft

50x

tg®

Porofix®

Mollelast® haft

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en Instructions for Use

Product description and performance characteristics

System for the treatment of venous leg ulcers over a period of approx. 12 weeks. For an effective compression bandage, apply 2 Rosidal® soft and 1 Rosidal® K (8 cm and 10 cm) each time. Replacement bandages are included in this set to allow washing after each dressing change. Use the consumables (tg®, Porofix® and Mollelast® haft) as needed (see bandaging instructions, items 1, 9, 16 and 17). If these are used up more quickly than expected (e.g. in cases of extensive oedema), the components can also be bought separately.

Product composition

- 1x tg® (tubular bandage), size 7, 20 m, 67% cotton, 33% viscose (Ref. 24005)
- 4x Rosidal® soft (foam bandage, white), 10 cm x 2 m x 0.2 cm, 100% polyurethane (Ref. 13016)
- 2x Rosidal® K (short stretch bandage, cream-coloured, 100% CO, 8 cm x 5 m, 100% cotton (Ref. 22201))
- 2x Rosidal® K (short stretch bandage, cream-coloured, 100% CO, 10 cm x 5 m, 100% cotton (Ref. 22202))
- 1x Porofix® (adhesive tape), 2.5 cm x 5 m, Support: 100% viscose, synthetic rubber adhesive (Ref. 34311)
- 1x Mollelast® haft (cohesive fixation bandage, white), 8 cm x 20 m, backing: 66% viscose, 34% polyamide, natural rubber adhesive (Ref. 30070)

Intended purpose

For venous leg ulcer compression therapy

Indications

Venous leg ulcers and their comorbidities, as well as chronic oedema (due to venous insufficiency)

Contraindications

Advanced peripheral arterial occlusive disease, decompensated heart failure, septic phlebitis, phlegmasia coerulea dolens, sensation disorders of the skin, known allergies and/or hypersensitivity to product components.

Side effects

In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions

- Incorrect bandaging can result in skin necrosis and pressure damage to nerves.
- Slight bluish tinge to the toes is normal; this should disappear as soon as the patient starts to walk.
- If the symptoms (pain and bluish tinge) do not disappear, the bandage must be removed and replaced with less compression.
- The patient must inform their doctor if sharply increasing pain occurs while wearing the bandage.
- With an Ankle Brachial Pressure Index (ABPI) greater than 0.8 or with an ABPI of 0.5–0.8 – only under medical supervision.
- In case of diabetes with microangiopathy, use only under medical supervision.

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Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 as follows: Dry cleaning waste code from Chapter 18 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

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