Suprasorb® CNP

Drainagepad
Dispositif de drainage
Drainage Pad
Almohadilla de drenaje
Almofada de drenagem
Pad di drenaggio
Drainagepad
Drænagepad
Dræningsdyna
Drøneringspad
Drenážna podložka
Dren z przystawką
Drenažna blazinica
Drenažni jastučić
Επίθεμα παροχέτευσης
Ταμρον de drenaj
Дренажна подложка
Drenažo tamponas
Drenaj ped
Drenážni polštářek
حشة التصريف

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1. Application:
a. Prior to use, visually inspect the sterile barrier system for breaches to the packaging integrity.
b. Open the peel pouch at the place provided for this purpose and remove the product, keeping it sterile.
c. Using sterile tweezers, hold the film that occludes the wound dressing (e. g., Suprasorb® F, Suprasorb® CNP EasyDress) as centrally as possible over the wound area and cut a round hole in the film about 3 cm in size with sterile scissors. (Figure 1).
d. To attach the Suprasorb® CNP Drainage Pad, pull the cover papers apart using the tabs and align the Suprasorb® CNP Drainage Pad so that the adhesive side faces the wound. (Figure 2).
e. Position the port of the Suprasorb® CNP Drainage Pad directly over the hole that you previously cut in the film. Keeping wrinkling to a minimum, press the adhesive film securely connected to the port on gently with your fingers so that the Suprasorb® CNP Drainage Pad adheres effectively and seals tightly. (Figure 3).
f. Then use the tabs to pull the cover paper off over the adhered film and tube of the Suprasorb® CNP Drainage Pad. (Figure 4).
g. Connect the Suprasorb® CNP Drainage Pad up to the Suprasorb® CNP P3 Therapy Unit using the connector. To do so, insert the exudate pouch tube of the Suprasorb® CNP P3 Therapy Unit two centimetres into the connector so that the end of the tube is visible in the window. (Figures 5 and 6).
h. Then lock the connector by pulling its ends apart. (Figures 7 and 8). The connections have now been securely established and you may begin the therapy.

2. Dressing change:
The Suprasorb® CNP Drainage Pad is changed at each dressing change.

3. Disposal:
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance - AVV). Recyclable packaging should be recycled according to the relevant national recycling system. Recyclable packaging, please observe any applicable national disposal regulations.

General instructions
In case of a serious incident, contact the manufacturer and the responsible health authorities.

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en Instructions for use

The instructions for use for the Suprasorb® CNP P3 system must also be complied with, in particular with regard to indications, contraindications, side effects, warnings and precautions.

Product description and features
The Suprasorb® CNP Drainage Pad consists of a flexible port with an exudate absorption zone in the centre and an integrated ventilation ring. Additional components include the adhesive film connected securely to the port, the four-lumen connecting tube (one lumen for exudate and three for ventilation), the tube clamp and the connector.

Product composition
Port: Polyvinyl chloride, cyanoacrylate adhesive
Film: Polyurethane, polyethylene, polyacrylate, paper coated with polydimethylsiloxane, polyester coated with polyacrylate
Connecting tube: Polyvinyl chloride
Tube clamp: Polypropylene
Connector: Polypropylene, acrylonitrile butadiene styrene and thermoplastic elastomers

Intended purpose
The Suprasorb® CNP Drainage Pad is placed over a wound dressing with the Suprasorb® CNP P3 Therapy Unit using the connector. To do so, insert the exudate pouch tube of the Suprasorb® CNP P3 Therapy Unit two centimetres into the connector so that the end of the tube is visible in the window. (Figures 5 and 6).
h. Then lock the connector by pulling its ends apart. (Figures 7 and 8). The connections have now been securely established and you may begin the therapy.

1. Aplicación:
a. Prior to use, visually inspect the sterile barrier system for breaches to the packaging integrity.
b. Open the peel pouch at the place provided for this purpose and remove the product, keeping it sterile.
c. Using sterile tweezers, hold the film that occludes the wound dressing (e. g., Suprasorb® F, Suprasorb® CNP EasyDress) as centrally as possible over the wound area and cut a round hole in the film about 3 cm in size with sterile scissors. (Figure 1).
d. To attach the Suprasorb® CNP Drainage Pad, pull the cover papers apart using the tabs and align the Suprasorb® CNP Drainage Pad so that the adhesive side faces the wound. (Figure 2).
e. Position the port of the Suprasorb® CNP Drainage Pad directly over the hole that you previously cut in the film. Keeping wrinkling to a minimum, press the adhesive film securely connected to the port on gently with your fingers so that the Suprasorb® CNP Drainage Pad adheres effectively and seals tightly. (Figure 3).
f. Then use the tabs to pull the cover paper off over the adhered film and tube of the Suprasorb® CNP Drainage Pad. (Figure 4).
g. Connect the Suprasorb® CNP Drainage Pad up to the Suprasorb® CNP P3 Therapy Unit using the connector. To do so, insert the exudate pouch tube of the Suprasorb® CNP P3 Therapy Unit two centimetres into the connector so that the end of the tube is visible in the window. (Figures 5 and 6).
h. Then lock the connector by pulling its ends apart. (Figures 7 and 8). The connections have now been securely established and you may begin the therapy.

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