
65 % PA, 35 % LYCRA®

70 % PA, 30 % LYCRA®

www.oeko-tex.com

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en Instructions for Use

Dear customer,

Thank you for choosing an L&R product.

Please read the instructions for use carefully. If you have any questions, please contact the person responsible for your treatment or your specialist shop.

Product description and performance characteristics

Actico® UlcerSys is a medical compression stocking system used to treat lower leg ulcers. The set includes an overstocking and two liners which also serve as donning aids to make it easier to put on and take off the overstocking.

The Actico® UlcerSys liner (refill pack) is part of the Actico® UlcerSys compression stocking system. However, the liner can also be used alone if only light compression is required.

If you are unable to don the compression stocking (overstocking) for physical reasons, a maximum of 2-3 liners can be worn over each other (only possible with the 10-14 mmHg liner, REF 32530-32539, 39740-39749).

The Actico® UlcerSys 16-18 mmHg liner (REF 454250-454259) is not meant to be worn in 2-3 layers over each other.

The refill pack (containing 3 liners) is available separately in either white or sand (only available

for the 10-14 mmHg version, REF 32530-32539, 39740-39749).
The set is available in two colour combinations – sand/white and black/sand.

Product composition

Overstocking: 65% PA polyamide, 35% Lycra®

Liner: 70% PA polyamide, 30% Lycra®

Intended purpose

The Actico® UlcerSys medical compression stocking system is intended for lower leg ulcer therapy over the treated wound and for the alleviation of venous symptoms.

Indications

- For venous leg ulcer therapy after oedema reduction
- Prevention of recurrence of chronic venous ulcers
- For stasis dermatitis (eczema, pigmentation, lipodermatosclerosis, atrophie blanche)
- For venous oedema
- For treatment of post-thrombotic syndrome
- For follow-up care for atrophy treatment and varicose vein surgery
- After vein removal for bypass surgery (cardiac surgery)

Contraindications

- Advanced peripheral arterial occlusive disease
- Congestive heart failure
- Septic phlebitis
- Phlegmasia cerulea dolens
- Diabetic ulcers
- Sensation disorders of the skin
- Known allergy and/or hypersensitivity to any of the product components

Side effects

Skin irritations and allergies can occur in very rare cases.

Warnings and precautions

Prior to application, cover any wound present with a wound dressing appropriate for the wound status.

- Due to the high resting pressure, remove the overstocking overnight.
- If sensation disorders of the extremities are present, consult the person responsible for your treatment prior to use.
- ABPI measurement: Actico® UlcerSys can be used with a value between 0.8 and 1.2. If the value is less than 0.8 or greater than 1.2, Actico® UlcerSys may only be used with extra caution and close monitoring, and only if recommended by the person responsible for treatment.

Usage information

- If severe oedema is present, the oedema must first be reduced using compression dressings (e.g. Rosidal® sys or Rosidal® TCS therapy systems).
- Once the lower leg ulcer has healed, it is essential to continue the compression treatment to prevent any recurrence. We recommend continuing the maintenance therapy either using the Actico® UlcerSys system or medical compression stockings specifically designed for this purpose.
- The Actico® UlcerSys should be applied by medically qualified skilled personnel the first time.
- The overstocking must be removed at night. The liner should be worn day and night.
- If toes become discoloured or unusual pain develops in the lower extremities when using the compression system, remove the stockings and immediately seek medical advice.
- The compression stockings must not come into contact with ointments, creams, oils or other similar substances.

Note regarding reuse

The product is intended exclusively for your use.

Donning instructions for the liner

- Place one hand in the liner and grab the heel of the liner with 2 fingers (Figure 1).
- Then turn the liner inside out to the point where the heel is being held (Figure 2).
- Pull the foot section over the foot until it is positioned over the forefoot (Figure 3).
- Pull the liner up around the heel. When doing so, ensure that the reinforced heel is correctly positioned on your heel (Figure 4).
- Now pull the liner over your ankle and calf up to below the knee (Figure 5).
- If you have drawn the liner up too high, simply push it back down to the ankle and then pull it up again using less force. Smooth out any wrinkles so that your stocking sits comfortably (Figure 6).

Donning instructions for the overstocking

- Place your hand in the overstocking and grab the heel of the overstocking with two fingers (Figure 7).
- Then turn the overstocking inside out to the point where the heel is being held (Figure 8).
- Pull the folded foot section over the foot until it is positioned over the forefoot (Figure 9).
- Now pull the overstocking up over the heel. Align the overstocking so that the inset stocking heel is positioned precisely over your heel. For optimal comfort in the toe area, pull the end of the liner forwards (Figure 10).
- Now pull the overstocking over your ankle and the calf up to below your knee and smooth out any wrinkles that form (Figure 11).
- Smooth out the overstocking evenly with both hands, using massage-like strokes until the stocking is free of any wrinkles to ensure that the pressure is evenly applied (Figure 12).

Washing instructions

Wash the liner using a standard mild detergent without fabric softener at 60°C (up to 50 machine washes) and the overstocking at 40°C (up to 50 machine washes). Lie flat to air dry. Do not wring out and do not tumble dry. Do not bleach, do not iron, do not dry clean. We recommend washing the liner after each dressing change. The liner can also be washed at 90°C in accordance with Austrian guidelines.

Warranty and storage instructions

The materials in compression stockings deteriorate over time and experience natural wear as a result of use. For this reason, the warranty period is 6 months under use.

Please store the compression stockings at room temperature and keep them away from direct sunlight, heat and moisture.

Disposal

In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Packaging labelled as recyclable should be taken into the respective national recycling systems.

General instructions

In case of a serious incident, contact the manufacturer and the responsible health authorities.

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