




Curafix® i.v. Control

Kanülenfixierpflaster mit Sichtfenster
Pansement de fixation de canule avec
fenêtre de contrôle
Cannula Retention Dressing with viewing
window
Apósito de fijación de cánulas con mirilla
Adesivo para fixação de cânulas com
janela de visualização
Cerotto per il fissaggio delle cannule con
finestra d'ispezione
Pleister met zichtvenster voor het fixeren
van canules

Kanylëfikseringsplaster med vindue
Kanylfixeringsplåster med transparent
inspektionsområde
Náplast k fixaci kanyl s okénkem
Fixačná náplast pre kanylu s okienkom
Ablakos kanülrögzítő tapasz
Plaster do mocowania kaniul z okienkiem
do kontroli miejsca wkłucia
Obliž za pritjevanje kanile s kontrolnim
okencem
Flaster za fiksiranje kanila s kontrolnim
prozorčićem
Gözlem pencereii kanül sabitleme bandı

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en Instructions for Use

Product description and features:
Curafix® i.v. Control is a retention dressing designed for peripheral access devices (in particular for catheters with injection ports). It consists of white non-woven material, adhesive material, film and cover paper. The dressing is provided with a separate padding cushion. In the top part of the product is a viewing window that allows monitoring of the cannula insertion site and the insertion channel. The dressing can be washed with water and iodine- or octenidine-containing disinfectants based on alcohol or water.

Product composition:
Polyester, polyacrylate, polyurethane, viscose, polypropylene, polyethylene, glassine paper with polydimethylsiloxane coating

Intended use:
For the retention of cannulas and catheters on the skin and as a mechanical barrier to protect the puncture site.

Indications:
Not applicable.

Contraindications:
Known allergy and/or hypersensitivity to any of the product components.

Side effects:
In very rare cases skin irritations and/or allergies may occur.

Warnings and precautions:
• The insertion site must be checked regularly for erythema, irritation and inflammation.
• If infection or other signs of inflammation occur (excessive warmth, oedema, erythema, pain), the decision to continue treatment must be made by healthcare professionals.
• Do not use if package is damaged or accidentally opened.

Instructions for use:
Application is performed by healthcare professionals.
1. Preparation:
After thorough cleaning and disinfection of the puncture site and the surrounding skin, the dressing is applied to the dried skin.
2. Use (applying the retention dressing – after puncture):
Before use, visually inspect the sterile barrier system for breaches to the integrity of the packaging.
a) Open the packaging at the point indicated and remove the product in a sterile manner (Fig. 1).
b) Before use, visually inspect the sterile barrier system for breaches to the integrity of the packaging.
c) Remove the large cover paper (Fig. 3).
d) Spread the dressing out at the wings and place the round incision around the injection port. The puncture site is now covered by the transparent viewing window (Fig. 4).
e) Now detach the 2 small cover papers from the wings and mould the dressing in place so that they slightly overlap (Fig. 5).
3. Dressing change:
The intervals at which the dressing must be changed are determined by the treating physician. These depend on the patient's individual skin condition and on the underlying medical findings.
Replace dressings which are soaked through, soiled or no longer adhere securely. Also replace the dressing if weeping pockets develop under the dressing.

Additional note:
To detach Curafix® i.v. Control, the use of dressing-loosening fluids (based on alcohol, saturated carbohydrites or siloxanes) is recommended.

4. Disposal:
In Europe, product waste can be assigned a waste code from Chapter 18 01 and packaging waste a waste code from Chapter 15 01 of the Directive on the European Waste Catalogue (Waste Catalogue Ordinance – AVV). Recyclable packaging should be recycled according to the relevant national recycling systems.

General instructions:
Curafix® i.v. Control is intended for single use and must not be resterilised. The product is sterile as long as the packaging stays unopened and undamaged.

In case of serious incident, contact the manufacturer and the responsible health authorities.

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