Gazin® RK

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en Instructions for Use

Product description
Gazin® gauze swabs RK, non-sterile, consist of 17- or 20-thread count cotton gauze compliant with EN 14079. The X-ray detectable thread contained in the swabs consists of barium sulphate filaments and is wound with polypropylene.

Product composition
100% cotton, X-ray detectable thread (blue, polypropylene, polyester thread, barium sulphate).

Intended purpose
Gazin® gauze swabs are intended to be used for wound cleansing or swabbing, for covering wounds as well as for absorbing body fluids during surgical procedures.

Usage information/Precautions
This product is non-sterile. Gazin® gauze swabs must be sterilised prior to surgically invasive use; resterilisation is not permitted. The product is intended for single use only.

Sterilisation instructions
1. Device requirements
Only sterilisers, which comply to the requirements of the norms EN 285 or EN 13060 and which use a pre-vacuum system for the sterilisation are allowed to be used. The use of the gravitational method is not allowed in the EU.
2. Sterilisation preparation
Before the product is sterilised, a suitable sterile barrier system must be selected, validated and used by the person responsible for the sterilisation (sterilising technician) in accordance with ISO 11607. The sterile barrier must be big enough for the product so that the seal is not under tension. The retail packaging is unsuitable for sterilisation.

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