Flat Knit THERAPY

Send completed form to Lohmann & Rauscher via:
Email: ordersaus@au.lrmed.com
Phone: 1300 572 869

VENOSAN®
Compression Stockings

VENOSAN® SOFT
VENOSAN® PRO
VENOSAN® STRONG


Lower Extremity Order Form

Date:

Therapist name:	Therapist phone:	Therapist email:
Patient name:	Patient phone:	Patient schema no:
Delivery details:		

Office use only:

Measurements




Hip/Waist

	Circumference (cm)	Length (cm)
cT		lGT
cH		lK1T
		lK2T

Foot/Leg

	Left circumference (cm)	Length (cm)	Right circumference (cm)
cG		lG	cG
cF		lF	cF
cE		lE	cE
cD		lD	cD
cC		lC	cC
cB1		lB1	cB1
cB		lB	cB
cY			cY
cA		lA	cA
		lZ	
		lA-medial	} must be filled in for oblique toe
		lA-lateral	

Garment types:



*short leg can finish at F or G

Additional garment styles are available. Contact L&R for details.

Product	Ccl.	Garment type	Qty. (pcs.)	Colour: beige	Colour: black
VENOSAN Soft	1				
VENOSAN Soft	2				
VENOSAN Pro	2				
VENOSAN Soft	3				
VENOSAN Pro	3				
VENOSAN Strong	3				
VENOSAN Strong	4				
VENOSAN Strong	4S				

Options

Bands
- ☐ 3 cm silicone band
- ☐ 5 cm silicone band
- ☐ 5 cm elastic band

Other
- ☐ open toe
- ☐ closed toe (please fill in lA and lZ)
- ☐ oblique toe (please fill in lA-medial and lA-lateral)
- ☐ T-heel
- ☐ horizontal fly
- ☐ open crotch
- ☐ full compression body part
- ☐ oblique G * *
- ☐ curved knee
- ☐ stocking in 2 parts (ie. AD + DG)
- ☐ stump finish
- ☐ waistband with button hole elastic

* * With AG and DG only possible in combination with 5 cm silicone or elastic band

Comments:

LR0295 0425

Flat Knit THERAPY

VENOSAN® Compression Stockings

VENOSAN® flat knit

Measuring Instructions - Lower Extremity

NOTE:

- All measurements in cm
- **c** = circumference, **l** = length
- Unless specified, all length measurements should be longitudinal and should not follow body contours

Leg Measurements

Measuring points B-D can be done with patient sitting or lying. Above E to be measured with patient standing.

cG Two fingers/4cm below Gluteal Fold for circumference

lG Two fingers/4cm below Gluteal Fold to ground for length

cF Mid Thigh for circumference

lF Mid Thigh to ground for length

cE Top of kneecap for circumference to avoid being tight at knee pit area

lE Middle of Kneecap to ground for length

cD Two fingers/4cm below Popliteal Fossa for circumference

lD Two fingers/4cm below Popliteal Fossa to ground for length

cC Widest part of Calf for circumference

lC Widest part of Calf to ground for length

cB1 Where Calf starts for circumference - approximately 6-8cm above point B

lB1 Where Calf starts to ground for length

cB Smallest part of the ankle, above Malleolus for circumference

lB Smallest part of the ankle, above Malleolus to ground for length

cY Around Instep and Heel with foot dorsiflexed for circumference

Foot Measurements

Around base of toes for circumference. Always in standing position, weight bearing.

Horizontal finish: Back of heel to base of toes for desired length

Oblique finish:

lA medial From the back of the heel to the base of large toe

lA lateral From the back of the heel to the base of small toe

Closed Toe: Back of heel to tip of longest toe (Total length of foot)

Panty Section Measurements

Ask patient to put their hands on their waistline

cT Waist circumference

cH Largest circumference around the hips or stomach

Measure the following length along contours of body

lGT Waist to 5cm/three fingers below Gluteal Fold

lK^{1}T Front length from centre of waist (T) to 5cm above Gluteal Fold applying gentle pressure

lK^{2}T Length at the back from the waist (T) over the buttocks to the Gluteal Fold while applying gentle pressure and following contours

Options

Silicone band: Helps keep garment in place.

Elastic band: makes top edge more firm. Does not prevent garment from slipping down. No length adjustments required for silicone or elastic bands.

T-Heel: can prevent rubbing or a crease indentation on the anterior ankle. Standard option in VENOSAN Strong garments.

Oblique G: Garment is finished on an angle.

Curved knee: Knee area knitted in a bent position.

Stocking in 2 parts: No additional measurements needed. The overlap will provide no additional pressure. Overlap is 7cm in length, if you do not state otherwise.

Stump finish: garment enclosed at end of garment to support base of stump.

Waist with button hole elastic: Waist band of panty section with button hole elastic to tighten around waist.

Please send the completed Order Form to Lohmann & Rauscher via:
Email: ordersaus@au.lrmed.com Phone: 1300 572 869

Scan to watch measuring video

www.Lohmann-Rauscher.com.au

LR0295 0425